


Elevate Integrated Actuator Module (EIAM) For Elevate Actuator Motor w/ MK6i Electronics



Elevate Integrated Actuator Module (EIAM) For Elevate Actuator Motor w/ MK6i Electronics

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1140036	Module, Elevate Integrated Actuator (EIAM) - MK6i (54" L)	1	1
2	A	1125685-NLA	NLA - Kit, MK5 SAC-E/ MK6i EIAM Mounting Hardware	1	1
3			Screw, Phillips Pan Head Tap (#10-32 x 1/2")	1	2
4			NLA - Bracket, Mounting, Black	1	1
5			NLA - Stud Plate, Quad Link	1	2

NOTE: A - Includes items 3-5