


Seat Frame Mounting Hardware (16"-22" Wide) Before 1/4/08 (CGE) - TDX SP/SR

Before 1/4/08 CGE was offered on narrow battery boxes only



Seat Frame Mounting Hardware (16"-22" Wide) Before 1/4/08 (CGE) - TDX SP/SR

Before 1/4/08 CGE was offered on narrow battery boxes only

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1144686	Kit, Seat Frame Mounting Hardware	1	1
2			NLA - Screw, Socket Head (5/16-18 x 3/4")	1	8
3			NLA - Washer, Ground (11/32 x 13/16 x 1/8")	1	8
4			Plate, Nut (5/16-18)	1	4
5		1142304	Screw, Hex Head SEMS w/ Patch (5/16-18 x 3/4")	1	8
6		1121721-NLA	NLA - Nut, Weld (5/16-18 x 1-1/8" Spacing)	1	4
7		1121710	NLA - Plate, Elevate Seat Mounting - Right	1	1
7		1121711-NLA	NLA - Plate, Elevate Seat Mounting - Left	1	1
NOTE: A - Includes items 2-4					