


# Seat Frame Mounting Hardware (16"-22" Wide) (CGTE/CGTRE) - TDX SP/SP-GT/SR



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1149787-NLA	NLA - Screw, Socket Flat Head (3/8-16 x 1-1/4")	1	4
2		1028769	Locknut (5/16-18)	1	4
3		1025708-NLA	NLA - Bolt, Hex Head Cap (1/4-20 x 3/4")	1	8
4		10502X001-NLA	NLA - Lockwasher (1/4")	1	8