



Battery Wire Harnesses For 22NF Batteries

For Narrow and Wide Battery Boxes (CGE/CGTE/CGTRE)



Battery Wire Harnesses For 22NF Batteries For Narrow and Wide Battery Boxes (CGE/CGTE/CGTRE)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1147093	Harness, Battery Wire For Elevate Option 22NF, Rear	1	1
2	A	1117223	Harness, Battery Wire For Elevate Option 22NF, Front	1	1
3		1096799-NLA	NLA - Screw, Hex Flange Whiz Lock (1/4-20 x 7/8")	1	2
4		16721X000	Locknut (1/4-20)	1	4
5		1095502-NLA	NLA - Package, Tie Wrap (11-1/2" L)	10	1
6		1095589-NLA	NLA - Battery Strap Assembly	1	3

NOTE: A - Items 3-5 are recommended when a new Battery Wire Harness is ordered