


Fender Shroud Assembly

TDX SR

Walking Beam, Motor and
Wheel Assembly shown for
reference only



Fender Shroud Assembly TDX SR

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1153415	Kit, Fender Shroud Assembly TDX SR - Right	1	1
2			Shroud, Fender, TDX SR - Right	1	1
3			Screw, Phillips Pan Head Tap (#10-32 x 1/2")	1	2
4			Decal, Push/Drive	1	1
1	A	1153416	Kit, Fender Shroud Assembly TDX SR - Left	1	1
2			Shroud, Fender, TDX SR - Left	1	1
3			Screw, Phillips Pan Head Tap (#10-32 x 1/2")	1	2
4			Decal, Push/Drive	1	1
5	B	1134240	Kit, Fender Shroud Mounting Bracket and Attaching Hardware	1	1
6			NLA - Screw, Socket Head, Cap (#10-32 x 3/4")	1	2
7			NLA - Bracket, Fender Shroud Mounting, TDX SR	1	1
8			NLA - Fastener, Loop with Adhesive Back (1" W X 1" L)	1	1
9			Locknut (#10-32)	1	2

NOTE: A - Includes items 2-4
B - Includes items 6-9