






# Battery Charger and Cables



## Battery Charger and Cables

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	F	1053161-NLA	NLA - Charger, Battery, 8 Amp, 115vac, Off Board Dual Mode Round Plug, Lester	1	1
1	A	1123249	Charger, Battery, 8 Amp. 115vac (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
1	A	1123249E-NLA	NLA - WARR RPLMT - Charger, Battery, 8 Amp, 115VAC (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
2	E	1089627-NLA	NLA - Cord, AC Power, 115VAC, 2 Prong Off-Board (11' 6" L)	1	1
2		1121086	Cord, AC Power, 115VAC, 3 Prong On-Board (6' 6-1/2" L)	1	1
3		60102824	Kit, Cable, Battery Extension	1	1
4	B	1146573-NLA	NLA - Cable, Drive Lockout Sensor, MK6i	1	1
5	C	1144926	Cable, Drive Lockout Sensor, MK6i	1	1
6	D	1142231	Cable, GB TT Interface	1	1
<p>NOTE: A - CTE charger is standard on-chair. Used on Gel Cell or AGM (Absorbed Glass Mat) batteries. For wet cell batteries refer to the Lester charger part number 1053161.</p> <p>B - Drive Lockout Sensor Cable used for chairs manufactured with ADJRECL Seat Frame and 4-Pole Motors only</p> <p>C - Sensor Cable used with item 6 Interface Cable for chairs manufactured with ADJRECL Manual Recline Seat Frame, MK6i Electronics and GB TT Motors</p> <p>D - Interface Cable used with item 5 Drive Lockout Sensor Cable for chairs manufactured with ADJRECL Manual Recline Seat Frame, MK6i Electronics and GB TT Motors</p> <p>E - No longer available as a service part. Please use 1121086 as its replacement.</p> <p>F - Use Battery Charger part number 1123249 as a replacement.</p>					