


Battery Box Tie-Down Bracket Mounting Hardware For TRRO/TRBKTS Option

Battery Box Tie-Down Bracket is a non-serviceable part. Call Technical Services at 1-800-832-4707.



Battery Box Assembly is shown for reference only

Battery Box Tie-Down Bracket Mounting Hardware For TRRO/TRBKTS Option

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1142304	Screw, Hex Head SEMS w/ Patch (5/16-18 x 3/4")	1	6