



# Actuator Motor, Elevate and Hardware with TRRO Option



**Elevate Hardware with TRRO Option  
For TDX SPREE with ADJASBA Seat**



# Actuator Motor, Elevate and Hardware with TRRO Option

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1122586	Sensor, Proximity Switch (30 1/4" L)	1	1
2		1154121-NLA	NLA - Kit, Elevate Hardware TRRO	1	1
3			NLA - Screw, Hex Head w/ Patch (5/16-18 x 7/8")	1	4
4			NLA - Screw, Hex Head w/ Patch (5/16-18 x 1/2")	1	4
5			NLA - Screw, Socket Head w/ Patch (5/16-18 x 1-1/2")	1	2
6			Screw, Socket Head (5/16-18 x 2-1/4")	1	4
7			Washer, Flat (5/16 x 5/8 x 1/16")	1	8
8			Locknut (5/16-18)	1	4
9			Clevis	1	1
10			Nut, Weld (5/16-18 x 1-3/8" Spacing)	1	2
11			Screw, Socket Flat Head w/ Patch (1/4-20 x 5/8")	1	2
12			Isolator, Rubber	1	1
13			Shim (1/64")	1	3
14	C B		Plug, TRRO/TRBKTS Tube - Right	1	1
14			Plug, TRRO/TRBKTS Tube - Left	1	1
15		1154206-NLA	NLA - Bracket, TRRO Support - Right	1	1
16		1152731	Kit, TRRO Elevate Actuator Motor	1	1
1			Sensor, Proximity Switch (30 1/4" L)	1	1
17			Actuator Motor Assembly, Elevate	1	1
18			Bracket, Lower Support	1	1
19			NLA - Screw, Hex Head w/ Patch (1/4-20 x 1-1/2")	1	6
20			Bracket, TRRO Support - Left	1	1
21			NLA - Screw, Socket Head w/ Patch (1/4-20 x 3/4")	1	2
22			Washer (1/4 x 5/8 x 3/32")	1	2
23			Tie Wrap (5-5/8" L)	1	2
24			Screw, Phillips Pan Head (#4-40 x 7/16")	1	4
25			NLA - Washer, Internal Tooth Lock (#4-40)	1	2
26			Nut, Hex Keps (#4-40)	1	2
27			Bracket, Limit Switch, Black	1	1
28			Sensor, Actuator Switch	1	1
24		1074128	Screw, Phillips Pan Head (#4-40 x 7/16")	1	4
25		1008525-NLA	NLA - Washer, Internal Tooth Lock (#4-40)	1	2
26		P0870X000	Nut, Hex Keps (#4-40)	1	2
28		1024001	Sensor, Actuator Switch	1	1

NOTE: A - Includes items 3-14  
B - Includes items 1 and 17-28  
C - Serial number of chair must be verified to ensure that TRRO option was originally ordered.