


Mercury Free Switches and Hardware

TDX SC

As of 7/1/2009, switches are mercury free.



Partial view of a Formula CG
Tilt Seating is shown for
reference only

Mercury Free Switches and Hardware TDX SC

As of 7/1/2009, switches are mercury free.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1140099	Switch, Drive Lockout Mercury Free Ball (MK6i) (50" L)	1	2
2	A	1154221-NLA	NLA - Kit, Switch Mounting Hardware CGT	1	1
3			Locknut (5/16-18)	1	2
4			NLA - Bracket, Cable/Wire Mounting, Black	1	2
5			Screw, Shoulder Hex Socket Head w/ Patch (3/8" x 3/8", 5/16-18)	1	2
NOTE: A - Includes items 3-5					