





Power Elevating Legrest - Articulating (MPCTMT/MPCTMT-L) TDX SC (Before 01/01/09)

NOTE: A complete set of Power Elevating Legrest-Articulating must be ordered thru Invacare's Equipment Configurator. Request Customer Service to review their reference document "Parts Pal" first before ordering.



**Power Elevating Legrest - Articulating
(MPCTMT/MPCTMT-L)**



**Power Elevating Legrest - Articulating
(MPCTMT/IFP)**

Power Elevating Legrest - Articulating (MPCTMT/MPCTMT-L) TDX SC (Before 01/01/09)

NOTE: A complete set of Power Elevating Legrest-Articulating must be ordered thru Invacare's Equipment Configurator. Request Customer Service to review their reference document "Parts Pal" first before ordering.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1147045	Kit, Power Elevating Legrest-Articulating Mtg. Hdwr. (MPCTMT/MPCTMT-L)	1	1
2			Screw, Socket Head Cap w/ Patch (1/4-20 x 1-1/2")	1	4
3			NLA - Spacer (1/4 x 3/8 x 3/16")	1	4
4			Cap, End, Side Rail, Black	1	2
5			Block, Threaded	1	2
6		1148708	Harness, Power Elevating Legrest-Articulating (MPCTMT/MPCTMT-L)	1	1
7	D	1153867-NLA	NLA - Harness, Tilt Actuator/MPCTMT Power Legrest and Sensor	1	1
	B,B1,C	1142951	Rail, Seat Frame (For 12 - 19" D Conventional)	1	1
	B,B1,C	1142952	Rail, Seat Frame (For 19 - 22" D Conventional)	1	1
<p>NOTE: A - Includes items 2-5 B - Not Shown. When a customer is ordering Power Elevating Legrest for the first time and a Contoura Seating System is used, the shorter seat rails used for Conventional Backs must be ordered also. B1 - This will prevent the legrest from being positioned to far forward C - Not Shown D - The Tilt Actuator/MPCTMT Power Legrest and Sensor Harness is used for Formula CGT on TDX SC models only</p>					