



# Pivot/ Slide Tube and Footplate

## P93ST / PW93ST

**NOTES:** Heel Loops are not included. See footrest and legrest accessories.  
For complete pivot/slide tube and footplate assemblies, see front riggings reference charts.



# Pivot/ Slide Tube and Footplate

## P93ST / PW93ST

NOTES: Heel Loops are not included. See footrest and legrest accessories.

For complete pivot/slide tube and footplate assemblies, see front riggings reference charts.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1120382	Tube, Pivot/Slide, Curved End, Black - Long (9-7/8")	1	1
1		1057852	Tube, Pivot/Slide, Curved End, Steel, Black (9-7/8")	1	1
2		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2
3		1027097	Plug Button, Black (3/4")	1	1
4		7930323-NLA	NLA - Package, Locknut (1/4-20)	25	1
5		7930366	Package, Socket Head Screw (1/4-20 x 2-1/4")	25	1
6		1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
6		1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
6		1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1