

# Aerosol Compressor - Stratos Portable (IRC1730)



**Stratos Portable Plus  
Shown for  
Reference Only**

<p>1</p>  A black power adapter with a two-prong electrical plug and a coiled black power cord with a standard connector.	<p>2</p>  A black flexible aerosol delivery tube with a standard connector at one end and a specialized nozzle at the other.	
<p>3</p>  A clear plastic aerosol chamber with a white nozzle and a coiled clear plastic tubing.	<p>4</p>  A clear plastic bag containing several white, circular filters.	<p>5</p>  A green and black carrying bag with a shoulder strap. The front panel features the Invacare logo and the text 'STRATOS PORTABLE'.

## Aerosol Compressor - Stratos Portable (IRC1730)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1156517-NLA	NLA - Adaptor, AC/DC Switching Portable	1	1
2		1133561-NLA	NLA - Adaptor, Car - Portable	1	1
3		MS2700S-NLA	NLA - Kit, Nebulizer Jet Plus	1	1
4		1156513	Package, Filter Replacement	5	1
5		IRC1722	Case, Carrying	1	1
	A	1133562-NLA	NLA - Cover, Filter	1	1
	A	1150755-NLA	NLA - Owner's Manual - Portable	1	1
	A	1150754-NLA	NLA - Instruction Sheet, Portable Battery	1	1

NOTE: A - Not Shown.