


# MKIV Swivel Joystick Hardware



## MKIV Swivel Joystick Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1093063	Spacer, Joystick Tube, Black	1	1
2		7930327-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 1-1/8")	25	3
3		1049962-NLA	NLA - Bracket, Swivel Joystick, MKIV	1	1
4		1049963-NLA	NLA - Bracket, Swivel Joystick with Hardware	1	1
5		1035986-NLA	NLA - Joystick Swivel Tube with Hardware, Black	1	1
6		1020551-NLA	NLA - Washer, Flat (1/4 x 3/8 x 1/16")	1	1
7		1005844-NLA	NLA - Bushing	1	2
8		1002511-NLA	NLA - Spring, Hitch Bracket	1	1
9		1023670-NLA	NLA - Knob, Clamping, Grey	1	1
10		1035990-NLA	NLA - Tube, Swivel Joystick, Black	1	1
11		1027095	Plug Button, Black (7/8")	1	1