

## Miscellaneous Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
		1114911-NLA	NLA - Knob, Adjustment, Black (5/16-18 x 1-1/2")	1	1
		40176X010-NLA	NLA - Pump, Tire	1	1
		See Chart	#1 Airline Style Positioning Strap	1	1
		See Chart	#1 Auto Style Push Button Positioning Strap	1	1
		See Chart	#1 Hook and Loop Style Positioning Strap	1	1

### #1 Auto Style Push Button Positioning Strap

Chair Width	Part Number	Seat/Chest Positioning Strap w/ Auto Style Buckle Length (Eyelet to Eyelet)
16" - 19" Wide	1004310	48" L
20" Wide / 18" Wide Van Seat	1034625	61" L
20" Wide Van Seat	1034622-NLA See Note(s) A	69" L

NOTE: A - For replacement contact The Aftermarket Group (TAG).

### #1 Airline Style Positioning Strap

Chair Width	Part Number	Seat/Chest Positioning Strap w/ Airline Style Buckle Length (Eyelet to Eyelet)
20" Wide / 18" Wide Van Seat	1034752-NLA See Note(s) A	61" L
20" Wide Van Seat	1034753-NLA	70" L

NOTE: A - 1034752 Airline Buckle Style is No Longer Available. Use Auto Style part number 1034625 as a replacement.

### #1 Hook and Loop Style Positioning Strap

Chair Width	Part Number	Seat/Chest Positioning Strap w/ Hook and Loop Style Buckle Length (Eyelet to Eyelet)
16" - 19" Wide	1034744 See Note(s) A	45" L
20" Wide	1034745 See Note(s) A	61" L

NOTE: A - Not available for Van Seat  
Hook and Loop Style Positioning Straps for power chairs used for chest positioning only.