


Frame Hardware



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Page	See Footnote	1	1
2		47471X074	Spring, Arm Lock	1	1
3		1028594	Pin, Arm Lock, Black	1	1
4		1090768	Kit, Lock Handle w/Rivets (Pair)	2	2
5		1027095	Plug Button, Black (7/8")	1	1
6	B	See Page	See Footnote	1	1

NOTE: A - See pages titled "Frames with Hardware".
B - See page titled "Back Canes".

NOTE: Items 2 thru 5 now come with frames.