


Power Back Reclining Mechanism



Power Back Reclining Mechanism

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1027095	Plug Button, Black (7/8")	1	6
2		60284X005-NLA	NLA - Screw, Socket Head (1/2-20 x 1-1/4")	1	2
3		1031658	Bracket, Controller Mounting, Black - Adult (18" W)	1	2
3		1031659	Bracket, Controller Mounting, Black - Narrow Adult (16" W)	1	2
4		1028514	Screw, Hex Head (5/16-24 x 1-1/2")	1	2
5		50200X000-NLA	NLA - Locknut (5/16-24)	1	8
6		1028496-NLA	NLA - Nut, Acorn (5/16-24)	1	6
7		1031891	Screw, Hex Head (5/16-24 x 2")	1	4
8		1031664-NLA	NLA - Screw, Shoulder Hex Socket Head (3/8 x 1")	1	2
9		1031886-NLA	NLA - Screw, Shoulder Hex Socket Head (1/2 x 1-3/4", 3/8-16)	1	2
10		60282X011-NLA	NLA - Screw, Shoulder Hex Socket Head (3/8 x 1-1/4", 5/16-18)	1	4
11		1031673-NLA	NLA - Washer, Nylon (35/64 x 1 x 1/16")	1	14
12		1031665	Screw, Shoulder Hex Socket Head (3/8 x 1-1/2", 5/16-18)	1	1
13		1029268	NLA - Rod, Motor Mount, Adult (18" W)	1	1
13		1029270-NLA	NLA - Rod, Motor Mount, Narrow, Adult (16" W)	1	1
14		14190X000-NLA	0	1	1
15		1025704	Washer, Flat (5/16 x 5/8 x 1/16")	1	5
16		1028769	Locknut (5/16-18)	1	3
17		50000X346-NLA	NLA - Mount, Spacer, Flange	1	1
18		1024220-NLA	NLA - Actuator Motor Assembly, Drive	1	1
19		1029262-NLA	NLA - Extender, Leg	1	2
20		1026163-NLA	NLA - Nut, Hex Jam. (1/2-20)	1	2
21		1029266-NLA	NLA - Rod, Rear Actuator, Adult (18" W)	1	1
21		1029264-NLA	NLA - Rod, Rear Actuator, Narrow, Adult (16" W)	1	1
22		1031885-NLA	NLA - Plate, Pivot (Black)	1	2
23		1031892-NLA	NLA - Spacer	1	1
24		1012371	Washer, Nylon (3/8 x 7/8 x 1/16")	1	6
25		47183X046-NLA	NLA - Spacer (11/32 x 5/8 x 1/2")	1	4
26		16325X000-NLA	NLA - Washer (3/8 x 13/16 x 1/16")	1	2
27		1028486	Locknut, Low Profile (3/8-16)	1	2
28		1031026	Screw w/ Lug (5/16-24 x 1-3/8")	1	2
29		1029255-NLA	NLA - Clamp, Bushing	1	2
30		1013149-NLA	NLA - Label, Caution	1	1
31		47447X009-NLA	NLA - Screw, Phillips (#8 x 3/8")	1	6
32	B,C	1037481-NLA	NLA - Headrest Extension Tube Assembly, Black	1	2
32	B,D	1037481-NLA	NLA - Headrest Extension Tube Assembly, Black	1	2
33		1029130-NLA	NLA - Tube, Recliner Back, Black - Right	1	2
33		1029131-NLA	NLA - Tube, Recliner Back, Black - Left	1	2
34	A	See Page	See Footnote	--	--

NOTE: A - See Frames with Hardware page.

B - Includes items 1 and 31

C - Without Heat-Shrink Tape - After 9/92

D - With Heat Shrink Tape - Before 10/92