



SuperPole with SuperBar / SuperPole with SuperBar Heavy Duty


PARTS DIAGRAM: SUPERPOLE™ with SUPERBAR™




PARTS DIAGRAM: POLE / CEILING ANCHOR PLATE ASSEMBLY



PARTS DIAGRAM: SUPERBAR



PARTS DIAGRAM: POLE / EXTENSION SCREW / BASE PLATE ASSEMBLY



SuperPole with SuperBar / SuperPole with SuperBar Heavy Duty

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP403	Ceiling Anchor Plate Assembly	1	1
1		SP403H DU	Top Plate, Ultra	1	1
1		SP403HD	Ceiling Anchor Plate Assembly, Heavy Duty	1	1
2		SP-P	Pole Assembly	1	1
2		SP-PHD	Pole Assembly, Heavy Duty	1	1
3	A	GR-101	Grip, 20"	1	1
4		SP404	SuperBar w/o the Castellated Collar	1	1
5		GR-100	Grip, 16"	1	1
6		SP402	Base Plate Assembly	1	1
7		HW323	Hitch Pin	1	1
8		HW318	Washer, Nylon	1	1
9		HW302	Clevis Pin	1	1
10		SP102	Pivot Bushing	1	2
11		HW100	Screw	1	2
12		HW300	Allen Key, Short Arm (3/32")	1	1
13		SP170	Spring Ring	1	1
14		HW112	Set Screws	1	4
15		SP411	Castellated Collar	1	2
16		HW320	Allen Key, Long Arm (5/32")	1	1
17		SP401	Cover Sleeve, Plastic	1	1
18		SP105	Threaded Collar	1	1
18		SP105HD	Threaded Collar, Heavy Duty	1	1
19		SP172	Jam Nut	1	1
20		SP107	Extension Screw	1	1
21		SP145	Washer, Plastic	1	1
Not Shown		SP-BX	Parts Box, SuperPole	1	1
Not Shown		SPB-BX	Parts Box, SuperBar	1	1
Not Shown		SPB-BXHD	Parts Box, SuperPole, Heavy Duty	1	1
Not Shown		SP501	Parts Pack, Super Bar	1	1

NOTE: A - Not available on Stainless Steel pole.