


Smart Rail

PARTS DIAGRAM: SMART-RAIL



Smart Rail

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SR403	Support Frame	1	1
2		SR402	Rail with Grip	1	1
3		HW301	Hex Key, Short Arm (1/8")	1	1
4		SR401	Latch Block	1	1
5		SR112	Support Frame	1	1
6		GR107	Traction Grips	1	4
7		SR109	Pivot Bushing	1	2
8		HW312	Cane Tips	1	2
9		SR400	Strap, Smart - Rail / Assista Rail	1	1
Not Shown		SR500	Parts Pack, Smart Rail	1	1