


# Assista-Rail

PARTS DIAGRAM: ASSISTA-RAIL



# Assista-Rail

| Item Number | Note | Part Number | Description                        | Quantity Per |          |
|-------------|------|-------------|------------------------------------|--------------|----------|
|             |      |             |                                    | Package      | Assembly |
| 1           |      | GR108       | Grip, Assista Rail Handle          | 1            | 1        |
| 2           |      | AST106-C    | Cross Bar                          | 1            | 1        |
| 3           |      | AST400      | Assista-Rail                       | 1            | 1        |
| 4           |      | GR-111      | Grips, Traction                    | 1            | 4        |
| 5           |      | SR400       | Strap, Smart - Rail / Assista Rail | 1            | 1        |
| Not Shown   |      | AST500      | Parts Pack, Assista-Rail           | 1            | 1        |