

Spoke Rear Wheels (w/o Handrims) (Before 9/2/11)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube Insert	Bearing	Hub Spacer
20SP	20 x 1-3/8"	Pneumatic	1115585-NLA See Note(s) C	1017519	1017520	40117X021	1015756
20SP	20 x 1-3/8"	Pneumatic w/ Flat Free	1115589-NLA See Note(s) A, C	1017519	1023322	40117X021	1015756
22SP	22 x 1-3/8"	Pneumatic	1115656-NLA See Note(s) C	1018544	1018545	40117X021	1015756
22SP	22 x 1-3/8"	Pneumatic w/ Flat Free	1115658-NLA See Note(s) A, C	1018544	1023323	40117X021	1015756
24SP	24"	Continental	1116705-NLA See Note(s) C	1026139-NLA	N/A	40117X021	1015756
24SP	24"	Pneumatic	1115660 See Note(s) C	40176X060	40176X061	40117X021	1015756
24SP	24"	Pneumatic w/ Flat Free	1116707 See Note(s) A, C	40176X060	1010906	40117X021	1015756
24SP	24"	Primo	1115662-NLA See Note(s) C	1026138	1032240	40117X021	1015756
24SP	24"	KIK	1114112 See Note(s) C	1062852	N/A	40117X021	1015756
24HPS	24"	Continental	1116713-NLA See Note(s) C	1026139-NLA	N/A	40117X021	1015756
24HPS	24"	Pneumatic	1116709-NLA See Note(s) C, E	40176X060	40176X061	40117X021	1015756
24HPS	24"	Pneumatic (Pair)	1116710	40176X060	40176X061	40117X021	1015756
24HPS	24"	Pneumatic w/ Flat Free	1116715 See Note(s) A, C	40176X060	1010906	40117X021	1015756
24HPS	24"	Primo	1116711 See Note(s) C	1026138	1032240	40117X021	1015756
24HPS	24"	KIK	1114138 See Note(s) C	1062852	N/A	40117X021	1015756
SPIN24	24"	Primo (Single #)	1112338-NLA See Note(s) D, B	1026138	1105941	40117X021	N/A
SPIN24	24"	Primo (1 Pair)	1112353-NLA See Note(s) D, B	1026138	1105941	1104350	N/A
559TS	25"	Primo	1116717-NLA See Note(s) C	1044947	1029471	40117X021	1015756
SPIN25	25"	Primo (Single #)	1112348-NLA See Note(s) D, B	1044947	1104326	1104350	N/A
SPIN25	25"	Primo (1 Pair)	1112354-NLA See Note(s) D, B	1044947	1104326	1104350	N/A
26HPS	26"	Continental	1115653-NLA See Note(s) C	1026143-NLA	N/A	40117X021	1015756
26HPS	26"	Pneumatic	1115652-NLA See Note(s) F, C	1026142	1026140	40117X021	1015756
26HPS	26"	Pneumatic	1116725-NLA	N/A	N/A	N/A	N/A
26HPS	26"	Primo	1115654 See Note(s) C	1029470	1029471	40117X021	1015756

Spoke Rear Wheels (w/o Handrims) (Before 9/2/11)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube Insert	Bearing	Hub Spacer
SPIN26	26"	Primo (Single #)	1112343-NLA See Note(s) D, B	1029470	1104326	1104350	N/A
SPIN26	26"	Primo (1 Pair)	1112355-NLA See Note(s) D, B	1029470	1104326	1104350	N/A
26SP	26"	Continental	1116721-NLA See Note(s) C	1026143-NLA	N/A	40117X021	1015756
26SP	26"	Pneumatic	1116719-NLA See Note(s) C	1026142	1026140	40117X021	1015756
26SP	26"	Primo	1116723-NLA See Note(s) C	1029470	1029471	40117X021	1015756
700SP	27"	Umma Gumma	1116727-NLA See Note(s) C	1045836-NLA	1045837-NLA	40117X021	1015756

NOTE: A - Stem Hole Plugs (Part Number 1030936) are used with Flat Free and Urethane Tires.
 B - As of 7/15/02 Spinergy Wheels now have 18 spokes instead of 30. For chairs built Before 7/15/02, it is recommended when one wheel is requested that the wheels be replaced in pairs.
 C - As of 11/08/02 these wheels now have a Black finish. For chairs built Before 11/08/02, it is recommended when one wheel is requested that the wheels be replaced in pairs.
 D - As of 11/17/08 Spinergy Wheels changed from black spoke, black hub to a black spoke, silver hub appearance. For chairs built Before 11/17/08, it is recommended when one wheel is requested that the wheels be replaced in pairs.
 E - Wheels must be replaced in pairs. See Part number 1116710 for replacement.
 F - No longer available as an individual service part. Wheel must be replaced in pairs, use 1116725 as its replacement.

Spoke Rear Wheels (w/o Handrims) (Before 9/2/11)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	1030936	Plug, Valve Stem	1	1
	A	See Chart	Wheels without Handrims - Spoke	1	1

NOTE: A - Not Shown