

Spoke Rear Wheels (w/ Handrims) (Before 9/2/11)

Wheel Style	Tire Type	Anodized Aluminum	Grey Plastic Coated	Black Plastic Coated	8 Vertical Projection	8 Oblique Projection
20SP	Pneumatic	1115584-NLA See Note(s) B	1049448-NLA	1115423-NLA See Note(s) B	N/A	N/A
20SP	Pneumatic w/ Flat Free	1115586 See Note(s) B	1049450-NLA	1115424-NLA See Note(s) B	N/A	N/A
22SP	Pneumatic	1115588-NLA See Note(s) B	1049452-NLA	1115591-NLA See Note(s) B	1115425-NLA See Note(s) B	N/A
22SP	Pneumatic with Flat Free	1115590-NLA See Note(s) B	1049456-NLA	1115593 See Note(s) B	1115426-NLA See Note(s) B	N/A
24SP	Pneumatic	1115594-NLA See Note(s) B	1049460-NLA	1115599 See Note(s) B	1115427-NLA See Note(s) B	1115428-NLA See Note(s) B
24SP	Pneumatic with Flat Free	1115597 See Note(s) B	1049465-NLA	1115601 See Note(s) B	1115429-NLA See Note(s) B	1115430-NLA See Note(s) B
24SP	Primo	1115602-NLA See Note(s) B	1049470-NLA	1115603-NLA See Note(s) B	1115431-NLA See Note(s) B	1115432-NLA See Note(s) B
24SP	Continental	1115606-NLA See Note(s) B	1049475-NLA	1115607-NLA See Note(s) B	1115433-NLA See Note(s) B	1115434-NLA See Note(s) B
24SP	KIK	1115645 See Note(s) B	1091332-NLA	1115646 See Note(s) B	1115447-NLA See Note(s) B	1115448-NLA See Note(s) B
24HPS	Pneumatic	1115622 See Note(s) B	1049496-NLA	1115623-NLA See Note(s) B	1115439-NLA See Note(s) B	N/A
24HPS	Pneumatic with Flat Free	1115625-NLA See Note(s) B	1049500-NLA	1115626-NLA See Note(s) B	1115440-NLA See Note(s) B	N/A
24HPS	Primo	1115628 See Note(s) B	1049504-NLA	1115629 See Note(s) B	1115441-NLA See Note(s) B	N/A
24HPS	Continental	1115631-NLA See Note(s) B	1049508-NLA	1115632-NLA See Note(s) B	1115442-NLA See Note(s) B	N/A
24HPS	KIK	1115649 See Note(s) B	1091337-NLA	1115650 See Note(s) B	1115449-NLA See Note(s) B	N/A
SPIN24	Primo	1112339 See Note(s) A, C	N/A	N/A	N/A	N/A
559TS	Primo	1115610 See Note(s) B	1049480-NLA	1115611-NLA See Note(s) B	1115435-NLA See Note(s) B	N/A
SPIN25	Primo	1112349-NLA See Note(s) A, C	N/A	N/A	N/A	N/A
26SP	Pneumatic	1115613-NLA See Note(s) B	1049484-NLA	1115614-NLA See Note(s) B	1115436-NLA See Note(s) B	N/A
26SP	Primo	1115616-NLA See Note(s) B	1049488-NLA	1115617-NLA See Note(s) B	1115437-NLA See Note(s) B	N/A
26SP	Continental	1115619-NLA See Note(s) B	1049492-NLA	1115620-NLA See Note(s) B	1115438-NLA See Note(s) B	N/A
26HPS	Pneumatic	1115636-NLA See Note(s) B	1049555-NLA	1115637 See Note(s) B	1115444-NLA See Note(s) B	N/A
26HPS	Primo	1115639-NLA See Note(s) B	1049559-NLA	1115640 See Note(s) B	1115445-NLA See Note(s) B	N/A
SPIN26	Primo See Note(s) A	1112344-NLA See Note(s) C, A	N/A	N/A	N/A	N/A
700SP	Umma Gumma	1115634 See Note(s) B	1049553-NLA	N/A	N/A	N/A

Spoke Rear Wheels (w/ Handrims) (Before 9/2/11)

NOTE: A - As of 7/15/02 Spinergy Wheels now have 18 spokes instead of 30. For chairs built Before 7/15/02, it is recommended when one wheel is requested that the wheels be replaced in pairs.
B - As of 11/8/02 these wheels now have a Black finish. For chairs built Before 11/8/02, it is recommended when one wheel is requested that the wheels be replaced in pairs.
C - As of 11/17/08 Spinergy Wheels changed from black spoke, black hub to a black spoke, silver hub appearance. For chairs built Before 11/17/08, it is recommended when one wheel is requested that the wheels be replaced in pairs.

When ordering a Spoke Wheel with Handrim attached, the standard handrim mounting process uses the outside holes of the handrim mounting tabs. To use the inside mounting holes to reduce the overall width of the chair, you must order the wheel and handrim assembled as a special order.