

## Spoke Rear Wheels (with Handrims)

| Wheel Style | Tire Type                | Anodized Aluminum               | Grey Plastic Coated | Black Plastic Coated         | 8 Vertical Projection        | 8 Oblique Projection         |
|-------------|--------------------------|---------------------------------|---------------------|------------------------------|------------------------------|------------------------------|
| 22SP        | Pneumatic                | 1115588-NLA<br>See Note(s) B    | 1049452-NLA         | 1115591-NLA<br>See Note(s) B | 1115425-NLA<br>See Note(s) B | N/A                          |
| 22SP        | Pneumatic with Flat Free | 1115590-NLA<br>See Note(s) B    | 1049456-NLA         | 1115593<br>See Note(s) B     | 1115426-NLA<br>See Note(s) B | N/A                          |
| 24SP        | Pneumatic                | 1115594-NLA<br>See Note(s) B    | 1049460-NLA         | 1115599<br>See Note(s) B     | 1115427-NLA<br>See Note(s) B | 1115428-NLA<br>See Note(s) B |
| 24SP        | Pneumatic with Flat Free | 1115597<br>See Note(s) B        | 1049465-NLA         | 1115601<br>See Note(s) B     | 1115429-NLA<br>See Note(s) B | 1115430-NLA<br>See Note(s) B |
| 24SP        | Primo                    | 1115602-NLA<br>See Note(s) B    | 1049470-NLA         | 1115603-NLA<br>See Note(s) B | 1115431-NLA<br>See Note(s) B | 1115432-NLA<br>See Note(s) B |
| 24SP        | Continental              | 1115606-NLA<br>See Note(s) B    | 1049475-NLA         | 1115607-NLA<br>See Note(s) B | 1115433-NLA<br>See Note(s) B | 1115434-NLA<br>See Note(s) B |
| 24SP        | KIK®                     | 1115645<br>See Note(s) B        | 1091332-NLA         | 1115646<br>See Note(s) B     | 1115447-NLA<br>See Note(s) B | 1115448-NLA<br>See Note(s) B |
| 24HPS       | Pneumatic                | 1115622<br>See Note(s) B        | 1049496-NLA         | 1115623-NLA<br>See Note(s) B | 1115439-NLA<br>See Note(s) B | N/A                          |
| 24HPS       | Pneumatic with Flat Free | 1115625-NLA<br>See Note(s) B    | 1049500-NLA         | 1115626-NLA<br>See Note(s) B | 1115440-NLA<br>See Note(s) B | N/A                          |
| 24HPS       | Primo                    | 1115628<br>See Note(s) B        | 1049504-NLA         | 1115629<br>See Note(s) B     | 1115441-NLA<br>See Note(s) B | N/A                          |
| 24HPS       | Continental              | 1115631-NLA<br>See Note(s) B    | 1049508-NLA         | 1115632-NLA<br>See Note(s) B | 1115442-NLA<br>See Note(s) B | N/A                          |
| 24HPS       | KIK®                     | 1115649<br>See Note(s) B        | 1091337-NLA         | 1115650<br>See Note(s) B     | 1115449-NLA<br>See Note(s) B | N/A                          |
| SPIN24      | Primo                    | 1112339<br>See Note(s) A, C     | N/A                 | N/A                          | N/A                          | N/A                          |
| 559TS       | Primo                    | 1115610<br>See Note(s) B        | 1049480-NLA         | 1115611-NLA<br>See Note(s) B | 1115435-NLA<br>See Note(s) B | N/A                          |
| SPIN25      | Primo                    | 1112349-NLA<br>See Note(s) A, C | N/A                 | N/A                          | N/A                          | N/A                          |
| 26SP        | Pneumatic                | 1115613-NLA<br>See Note(s) B    | 1049484-NLA         | 1115614-NLA<br>See Note(s) B | 1115436-NLA<br>See Note(s) B | N/A                          |
| 26SP        | Primo                    | 1115616-NLA<br>See Note(s) B    | 1049488-NLA         | 1115617-NLA<br>See Note(s) B | 1115437-NLA<br>See Note(s) B | N/A                          |
| 26SP        | Continental              | 1115619-NLA<br>See Note(s) B    | 1049492-NLA         | 1115620-NLA<br>See Note(s) B | 1115438-NLA<br>See Note(s) B | N/A                          |
| 26HPS       | Pneumatic                | 1115636-NLA<br>See Note(s) B    | 1049555-NLA         | 1115637<br>See Note(s) B     | 1115444-NLA<br>See Note(s) B | N/A                          |
| 26HPS       | Primo                    | 1115639-NLA<br>See Note(s) B    | 1049559-NLA         | 1115640<br>See Note(s) B     | 1115445-NLA<br>See Note(s) B | N/A                          |
| SPIN26      | Primo<br>See Note(s) A   | 1112344-NLA<br>See Note(s) C, A | N/A                 | N/A                          | N/A                          | N/A                          |
| 700SP       | Umma Gumma               | 1115634<br>See Note(s) B        | 1049553-NLA         | N/A                          | N/A                          | N/A                          |

## Spoke Rear Wheels (with Handrims)

NOTE: A - As of 7/15/02 Spinergy Wheels now have 18 spokes instead of 30. For chairs built Before 7/15/02, it is recommended when one wheel is requested that the wheels be replaced in pairs.  
B - As of 11/8/02 these wheels now have a Black finish. For chairs built Before 11/8/02, it is recommended when one wheel is requested that the wheels be replaced in pairs.  
C - As of 11/17/08 Spinergy Wheels changed from black spoke, black hub to a black spoke, silver hub appearance. For chairs built Before 11/17/08, it is recommended when one wheel is requested that the wheels be replaced in pairs.

When ordering a Spoke Wheel with Handrim attached, the standard handrim mounting process uses the outside holes of the handrim mounting tabs. To use the inside mounting holes to reduce the overall width of the chair, you must order the wheel and handrim assembled as a special order.