


MK5 Display



| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|------|-------------|---|--------------|----------|
| | | | | Package | Assembly |
| 1 | | 1118759-NLA | NLA - Display, MK5 with Switches | 1 | 1 |
| 2 | | 1093071 | Clamp, Fixed Height Joystick Tube | 1 | 1 |
| 3 | A | 1139422 | Clamp, Height Adjustable Joystick Tube (ARM250) | 1 | 1 |
| 4 | | 1049594-NLA | NLA - Tube, Flex with Screws | 1 | 1 |
| 5 | | 1052041 | Screw, Phillips Pan Head (#4-40 x 3/8") | 1 | 4 |

NOTE: A - Introduced as an option (ARM250) and can be used to replace a Fixed Height Joystick Clamp