


Multi-Function Frame

Motor Mounting Hardware

(When Using Power Reclining Back)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1014370-NLA	NLA - Screw, Socket Head (1/2-20 x 3-1/2")	1	1
2		1031673-NLA	NLA - Washer, Nylon (35/64 x 1 x 1/16")	1	1
3		1013930-NLA	NLA - Spacer (1/2 x 3/4 x 1-3/4")	1	1
4		90070X000-NLA	NLA - Nut (1/2-13)	1	1
5		1014371	Lockwasher (1/2")	1	1
6		7930323-NLA	NLA - Package, Locknut (1/4-20)	25	1
7		1013925-NLA	NLA - Spacer (1/2 x 3/4 x 7/16")	1	1
8		1013923-NLA	NLA - Mounting Plate	1	1
9		60274X011-NLA	NLA - Screw, Hex Cap (1/4-20 x 1/2")	1	1
10		60282X006-NLA	NLA - Bolt, Socket Head Shoulder (1/2 x 1-3/4")	1	1
11		1000134-NLA	NLA - Screw, Socket Button Head (1/2-13 x 1")	1	1
12		48323X013	Washer (5/8 x 7/8 x 1/8")	1	1