


# Back Canes

(Contoura Back, ABS & Metal Curved Back and Tilt w/ BPO Option)



ADJASBA Seat Frame and Seat Pan  
shown for reference only

# Back Canes

(Contoura Back, ABS & Metal Curved Back and Tilt w/ BPO Option)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Back Cane Assemblies	1	1
2		1135506	Package, Foam Hand Grip (9" L)	10	1
3		1134426	Package, Plug Button, Black (3/4")	10	1

NOTE: A - Includes the Back Cane and items 2 and 3

## #1 Back Cane Assemblies

Back Height	Back Cane-Push Handle Straight Right	Back Cane-Push Handle Straight Left
16"	1133793	1133794
17"	1133795	1133796
18"	1143334 See Note(s) C	1143335 See Note(s) C
19"	1133799	1133800
20" See Note(s) C	1133801 See Note(s) D, B	1133802 See Note(s) D, B
21"	1136011 See Note(s) A	1136012 See Note(s) A
22"	1136013 See Note(s) A	1136014 See Note(s) A
23"	1136015 See Note(s) A	1136016 See Note(s) A
24"	1136017 See Note(s) A	1136018 See Note(s) A

NOTE: A - 21"-24" high straight back canes with push handles are only offered On-Chair when a BPO Option is ordered  
 B - All (20"-24" high ABS Curved Backs) or (20"-26" high Metal Curved Backs) ordered On-Chair will use the 20" high straight cane w/ push handle  
 C - Wheelchairs manufactured before 5/15/09 with a Contoura Back 25" or 29" high and wheelchairs manufactured after 5/14/09 with a Curve Contoura Back 22" high will use the 18" high Straight Back Cane w/ Push Handle.  
 D - Wheelchairs manufactured after 5/14/09 with a Curve Contoura Back 25" or 28" high will use the 20" high Straight Cane w/ Push Handle.

NOTE: Back Cane Assemblies include Foam Grips and Plug Buttons where shown