




Seat Positioning Straps



Seat Positioning Straps

| Item Number | Note | Part Number | Description | Quantity Per | |
|---|------|-------------|---|--------------|----------|
| | | | | Package | Assembly |
| 1 | A | 1004310 | Strap, Positioning, Seat/Chest w/ Auto Style Push Button Buckle (Eyelet to Eyelet Length 48") | 1 | 1 |
| 1 | B | 1034625 | Strap, Positioning, Seat/Chest w/ Auto Style Push Button Buckle (Eyelet to Eyelet Length 61") | 1 | 1 |
| 1 | B,E | 1034752-NLA | NLA - Strap, Positioning, Seat/Chest w/ Airline Style Buckle (Eyelet to Eyelet Length 61") | 1 | 1 |
| 1 | C | 1034745 | Strap, Positioning, Seat/Chest w/ Hook and Loop Style Buckle (Eyelet to Eyelet Length 61") | 1 | 1 |
| 2 | | 1037610 | Screw, Hex Head (1/4-20 x 1-5/8") | 1 | 2 |
| 3 | | 1025736 | Locknut (1/4-20) | 1 | 2 |
| 4 | D | 1026298 | Package, Screw with Washer (#10 x 7/8") (Pair) | 2 | 1 |
| NOTE: A - 16"-19" Wide B - 20"-22" Wide C - 16"-22" Wide D - Hardware used for mounting the Chest Positioning Strap E - 1034752 Airline Buckle Style is No Longer Available. Use Auto Style part number 1034625 as a replacement. | | | | | |