


Base Assembly Components


Model 1200, 1300, 1500, 1600, 1800 Service Parts

Standard, Narrow and C-Base Lift Assemblies

1




2




3



6




8



5

5

4



7

7

9



SICC_033A

Base Assembly Components

Model 1200, 1300, 1500, 1600, 1800 Service Parts

Standard, Narrow and C-Base Lift Assemblies

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		TRP1000A-NLA	NLA - Standard Base, Almond	1	1
1		TRP1000W	Standard Base, White	1	1
2		TRP1000CA-NLA	NLA - C-Base, Almond	1	1
2		TRP1000CW-NLA	NLA - C-Base, White	1	1
3		IHTRP1029-NLA	NLA - Caster, Brake - Base	1	1
4		IHTRP1028-NLA	NLA - Caster, Steering - Base	1	1
5		IHTRP1062B	Bolt, Caster (M8 x 25mm)	1	2
6		IHTRP1027-NLA	NLA - Caster, Front - Base	1	2
7		IHTRP1063-NLA	NLA - Bolt, Stainless, Front Caster (7/16-18 x 1")	1	2
8		IHTRP1001A-NLA	NLA - Motor Box, Almond	1	1
8		IHTRP1001W-NLA	NLA - Motor Box, White	1	1
9		IHTRP1065-NLA	NLA - Bolt, Motor Box, Stainless (1/4-20 x 3/8")	1	4