


# Model 1800 - C Base Seat Assembly



SICC\_041B

## Model 1800 - C Base Seat Assembly

| Item Number                                                                                                                 | Note | Part Number     | Description                                 | Quantity Per |          |
|-----------------------------------------------------------------------------------------------------------------------------|------|-----------------|---------------------------------------------|--------------|----------|
|                                                                                                                             |      |                 |                                             | Package      | Assembly |
| 1                                                                                                                           | A,B  | IHTRV1801       | Model 1800 Seat Assembly                    | 1            | 1        |
| 2                                                                                                                           | C    | TRP1802A-NLA    | NLA - High Back Commode Seat - Almond       | 1            | 1        |
| 2                                                                                                                           | C    | TRP1802W-NLA    | NLA - High Back Commode Seat - White        | 1            | 1        |
| 3                                                                                                                           | C    | TRV1803A-NLA    | NLA - High Back Solid Seat - Almond         | 1            | 1        |
| 3                                                                                                                           | C    | IHTRP1803W-NLA  | NLA - High Back Solid Seat - White          | 1            | 1        |
| 4                                                                                                                           | C    | IHTRP1802PA-NLA | NLA - High Back Solid Seat w/ Pail - Almond | 1            | 1        |
| 4                                                                                                                           | C    | IHTRP1802PW-NLA | NLA - High Back Solid Seat w/ Pail - White  | 1            | 1        |
| 5                                                                                                                           |      | IHTRP1804-NLA   | NLA - Commode Pail Package                  | 1            | 1        |
| 6                                                                                                                           |      | IHTRP9702-NLA   | NLA - Footman Loop                          | 1            | 1        |
| 7                                                                                                                           |      | IHPAR9721C      | Harness Strap                               | 1            | 1        |
| 8                                                                                                                           |      | IHTRP1302C-NLA  | NLA - Seat Hanger Rod                       | 1            | 1        |
| 9                                                                                                                           |      | IHTRP1806-NLA   | NLA - C Clamp - 1" Tube                     | 1            | 1        |
| 10                                                                                                                          |      | IHTRP1034C      | Tube, Foam (7/8" x 8")                      | 1            | 1        |
| 11                                                                                                                          |      | IHTRP1807-NLA   | NLA - Locknut (1/4-20)                      | 1            | 1        |
| 12                                                                                                                          |      | IHTRP1034A-NLA  | NLA - Tube, Foam (7/8" x 24")               | 1            | 2        |
| 13                                                                                                                          |      | IHTRP1304LC-NLA | NLA - Left Arm Rest, Seat Assembly          | 1            | 1        |
| 14                                                                                                                          |      | IHTRP1304RC-NLA | NLA - Right Arm Rest, Seat Assembly         | 1            | 1        |
| 15                                                                                                                          |      | IHTRP1037-NLA   | NLA - End Cap (7/8")                        | 1            | 2        |
| 16                                                                                                                          |      | IHTRP1311       | Roll Pin (3/16 x 3/4")                      | 1            | 2        |
| 17                                                                                                                          |      | IHTRP1310-NLA   | NLA - Roll Pin (3/16 x 1")                  | 1            | 2        |
| NOTE: A - Specify color and solid or commode when ordering.<br>B - Includes Items 2 thru 17.<br>C - Includes Items 6 and 7. |      |                 |                                             |              |          |