

Footrest and Legrest Accessories

| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|------|-------------|--|--------------|----------|
| | | | | Package | Assembly |
| | A,B | 1035902 | Strap, Calf w/ Hook and Loop (31" L) | 1 | 1 |
| | A,C | 1032243 | Strap, Calf w/ Hook and Loop (36" L) | 1 | 1 |
| | A,D | 1026987 | Guard, Impact (8" L) | 1 | 2 |
| | A | See Chart | Heel Loop Attaching Hardware (1650, 1450, & 1350 Footplates) | 1 | 1 |
| | A | See Chart | Heel Loops and Hardware (1650, 1450, & 1350 Footplates) | 1 | 1 |
| | A | See Chart | Heel Loops and Hardware (AT5543 Footplates) | 1 | 1 |

NOTE: A - Not Shown
 B - 18-22" Wide chairs
 C - 23-24" Wide chairs
 D - For 70 deg. Tapered Footrest.

Heel Loops and Hardware

For AT5543 Footplates

| Footplate Type | Part Number | Part Number Consists of |
|-------------------|---|---|
| AT5543 Adj. Angle | 1160723 Kit, Heel Loop w/ Hardware, Adjustable Angle Footplate (Pair) | (2) Heel Loops (4) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (4) Spacer (1/4 x 5/16 x 1-3/32") (4) Locknut (1/4-20) |

Heel Loops and Hardware

For 1650, 1450, & 1350 Footplates

| Footplate Type | Chair Width | Heel Loop Single # | Heel Loop Pair # | Attaching Hardware Pair # |
|--|-------------|--------------------|------------------|---------------------------|
| Standard Aluminum (1650) Deluxe Aluminum (1450) | 18-24" Wide | 1017621 | 1035917-NLA | 1035907 |
| X-Large (1350) | 18-24" Wide | 1027000 | 1035929 | 1055241 |

Heel Loop Attaching Hardware

For 1650, 1450, & 1350 Footplates

| Footplate Type | Part Number | Part Number Consists of |
|----------------|-------------|-------------------------|
|----------------|-------------|-------------------------|

Footrest and Legrest Accessories

Heel Loop Attaching Hardware

For 1650, 1450, & 1350 Footplates

| Footplate Type | Part Number | Part Number Consists of |
|--|-------------|--|
| Standard Aluminum (1650) Deluxe Aluminum (1450) | 1035907 | 2 each Screw, PPH (1/4-20 x 2-1/2") 2 each Spacer (1/4 x 11/32 x 1-3/4") 2 each Locknut (1/4-20) |
| X-Large (1350) | 1055241 | 2 each Screw, PPH (1/4-20 x 2-1/4") 2 each Spacer (1/4 x 11/32 x 1-3/4") 2 each Locknut (1/4-20) |