

## Footrest and Legrest Accessories

**NOTE: Heel loops for Power Articulating Legrest located on Power Articulating Legrest page.**

### Toe Loops (1345)

Footplate Type	Toe Loop w/o Hardware Single #	Toe Loop w/o Hardware Pair #	Attaching Hardware Part Number	Attaching Hardware Consists of
1651, 1650, 1450, 1350, Footplates	04572X025	1037087	1037089	(4) Screw, Phillips Pan Head (#10-32 x 7/8"), (4) Washer (1/4 x 3/4 x 3/32") (4) Locknut (#10-32)
AT5543 Footplates	04572X025	1037087	1104749	(4) Screw, Phillips Pan Head (#10-32 x 1/2") (4) Washer (1/4 x 3/4 x 3/32")

### Heel Loops (1346)

Footplate Type	Heel Loop Single #	Heel Loop Pair #	Heel Loop with Ankle Strap - Right Single #	Heel Loop with Ankle Strap - Left Single #	Heel Loop with Ankle Strap Pair #	Attaching Hardware Pair #
Composite (1651) 16" Wide Tapered	1028412-NLA	1035912-NLA	1035908-NLA	1035909-NLA	1035906	1035931
Composite (1651) 18" Wide Tapered	1017621	1035917-NLA	1035919	1035920	1035918	1035931
Composite (1651) 20" Wide Tapered	1017621	1035917-NLA	1035932-NLA	1035933-NLA	1035930	1035931
Standard Aluminum (1650) 14" Wide	1025234-NLA	1035905	1035908-NLA	1035909-NLA	1035906	1035907
Standard Aluminum (1650) 16-20" Wide	1017621	1035917-NLA	1035919	1035920	1035918	1035907
Deluxe Aluminum (1450) 16-20" Wide	1017621	1035917-NLA	1035919	1035920	1035918	1086112-NLA

## Footrest and Legrest Accessories

**NOTE: Heel loops for Power Articulating Legrest located on Power Articulating Legrest page.**

### Heel Loops (1346)

Footplate Type	Heel Loop Single #	Heel Loop Pair #	Heel Loop with Ankle Strap - Right Single #	Heel Loop with Ankle Strap - Left Single #	Heel Loop with Ankle Strap Pair #	Attaching Hardware Pair #
X-Large Aluminum (1350) 14" Wide	1025234-NLA	1035905	1035908-NLA	1035909-NLA	1035906	1055241
X-Large Aluminum (1350) 16-20" Wide	1017621	1035917-NLA	1035919	1035920	1035918	1055241
Adjustable Angle (AT5543)	1006327-NLA See Note(s) B	1160723 See Note(s) A	1038229	1038229	1160724 See Note(s) A	N/A
<b>NOTE:</b> A - These Assemblies include the Attaching Hardware. B - No longer available. Use part number 1160723 as replacement.						

## Footrest and Legrest Accessories

**NOTE: Heel loops for Power Articulating Legrest located on Power Articulating Legrest page.**

### Heel Loop Hardware

Footplate Type	Part Number	Consists of:
Adjustable Angle (AT5543)	1160723	(2) Heel Loops (4) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (4) Spacer (4) Locknut (1/4-20)
Adjustable Angle (AT5543)	1160724	(2) Heel Loops w/ Ankle Straps (4) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (4) Spacer (4) Locknut (1/4-20)
Composite 16-20" W (1651)	1035931	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Hex Nut w/ Nylon Cap (1/4-20) (2) Flat Washer (1/4 x 5/8 x 1/16")
X-Large Aluminum 14-20" W (1350)	1055241	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Locknut (1/4-20")
Standard Aluminum 14-20" W (1650)	1035907	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/2") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Locknut (1/4-20")
Deluxe Aluminum 16-20" W (1450)	1086112-NLA	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/2") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Locknut (1/4-20")

### Miscellaneous Accessories

Part Number	Description
1026987 See Note(s) B	Impact Guard (8") (IMPACTG)
1143050	Non-Removable Device (1820)
1032003	"H" Heel Strap (9342) 14" Wide
1032002	"H" Heel Strap (9342) 16-20" Wide
1035902	Calf Strap (1337)

# Footrest and Legrest Accessories

NOTE: Heel loops for Power Articulating Legrest located on Power Articulating Legrest page.

## Miscellaneous Accessories

Part Number	Description
8881006725-NLA  See Note(s) A	Footrest Support Panel (1348) 14-16" Wide
8881006726  See Note(s) A	Footrest Support Panel (1348) 18" Wide
8881006727-NLA  See Note(s) A	Footrest Support Panel (1348) 20" Wide
8881035816-NLA  See Note(s) A	Legrest Support Panel (1348) 14-16" Wide
8881035817-NLA  See Note(s) A	Legrest Support Panel (1348) 18" Wide
8881035818-NLA  See Note(s) A	Legrest Support Panel (1348) 20" Wide
NOTE: A - Specify upholstery color when ordering. See page titled "Frame Finishes, Touch-Up Paints and Upholstery Colors." B - For 70 deg. Tapered Footrest.	
NOTE: These items are not shown	