


# AT5543 Adjustable Angle Footplates

## PW93 and PW93ST Swingaway Supports

NOTE: Heel loops are not included. See footrest and legrest accessories.



## AT5543 Adjustable Angle Footplates PW93 and PW93ST Swingaway Supports

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	PW93 Pivot/Slide Tube and Footplate Assemblies	1	1
1	A	See Chart	PW93ST Pivot/Slide Tube and Footplate Assemblies	1	1
2	B	1113573	Footplate, Adjustable Angle, Large, Black (5-1/4" W x 8" L)	1	1
2	B	1113574	Footplate, Adjustable Angle, Extra Large, Black (6-1/4" W x 8" L)	1	1
3		1058517	Clamp, Footplate	1	1
4		1057855	Tube, Pivot and Slide, Straight, Steel, Black (9-7/8" L)	1	1
4		1057856	Tube, Pivot/Slide, Straight, Steel, Black (11-7/8" L)	1	1
5		1048816	Screw, Nylon Adjustment (1/2-13 x 7/8")	1	1
6		1052102	Hinge, Footplate - Right	1	1
6		1052103	Hinge, Footplate - Left	1	1
7		1029545	Screw, Socket Flat Head (1/4-20 x 3/4")	1	2
8		1082686	Nut, Barrel, Hex Head (1/4-20 x 11/16")	1	2
9		7930321-NLA	NLA Package, Screw, Socket Button Head (1/4-20 x 1-7/8")	25	1
10		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	1
11		1026900	Washer (1/4 x 23/32 x 1/16")	1	1
12		7930313	Package, Washer, Nylon (5/16 x 7/8 x 1/8")	25	1
13		1026907	Washer (1/4 x 1/2 x 1/16")	1	1
14		7930333	Locknut (1/4-20)	25	1

NOTE: A - Includes Items 2 thru 14.  
B - This footplate was revised 9/2/2008 with rounder corners and accessory holes changed from 8 to 4. Recommended to replace footplates in pairs to match sides.

### Item #14 PW93 Pivot/Slide Tube and Footplate Assemblies

Chair Width	16" and Wider	16" and Wider
Footplate Type	6-1/4" X-Large Right	6-1/4" X-Large Left
PW93 Swingaway Supports 9-7/8" Slide Tube	1062756	1062757
PW93 Swingaway Supports 11-7/8" Slide Tube	1062760	1062761

### Item #14 PW93ST Pivot/Slide Tube and Footplate Assemblies

Chair Width	16" and 17" Wide	16" and 17" Wide	18" and Wider	18" and Wider
Footplate Type	5-1/4" Large Right	5-1/4" Large Left	6-1/4" X-Large Right	6-1/4" X-Large Left
PW93ST Swingaway Supports	1060444	1060445	1062756	1062757

# AT5543 Adjustable Angle Footplates PW93 and PW93ST Swingaway Supports

## Item #14 PW93ST Pivot/Slide Tube and Footplate Assemblies

Chair Width	16" and 17" Wide	16" and 17" Wide	18" and Wider	18" and Wider
PW93ST Swingaway Supports with 2" Longer Slide Tube	1062758	1062759	1062760	1062761