

Footrest and Legrest Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	1032243	Strap, Calf w/ Hook and Loop (36" L)	1	1
	A,B	1032215	H-Strap (26" L)	1	1
	A,C	1032216	H-Strap (31" L)	1	1

NOTE: A - Not Shown
 B - 16" wide seats
 C - 17" wide seats and wider

Footrest and Legrest Accessories (Not Including AT5543)

Footplate Type	Chair Width	Heel Loop Single #	Heel Loop Pair #	Heel Loop w/ Ankle Strap Single # Rt.	Heel Loop w/ Ankle Strap Single # Lt.	Heel Loop w/ Ankle Strap Pair #	Attaching Hardware Pair #
Composite Small (1651)	16-17" W Tapered	1028412-NLA	1035912-NLA	1035908-NLA	1035909-NLA	1035906	1035931
Composite Medium (1651)	16" W or 18-19" W Tapered	1017621	1035917-NLA	1035919	1035920	1035918	1035931
Composite Large (1651)	17-24" W or 20-24" W Tapered	1017621	1035917-NLA	1035932-NLA	1035933-NLA	1035930	1035931
Standard Aluminum (1650)	16-24" W	1017621	1035917-NLA	1035919	1035920	1035918	1035907
Deluxe Aluminum (1450)	16-24" W	1017621	1035917-NLA	1035919	1035920	1035918	1086112-NLA
X-Large Aluminum (1350)	16-24" W	1017621	1035917-NLA	1035919	1035920	1035918	1055241

Attaching Hardware (Not Including AT5543)

Attaching Hardware Part Numbers	Attaching Hardware Part Numbers
1035931	1055241
(2) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Hex Nut w/ Nylon Cap (1/4-20) (2) Flat Washer (1/4 x 5/8 x 1/16")	2 ea. Screw, Phillips Pan Head (1/4-20x2-1/4") 2 ea. Spacer (1/4x11/32x1-3/4") 2 ea. Locknut (1/4-20)

Footrest and Legrest Accessories

Attaching Hardware (Not Including AT5543)

Attaching Hardware Part Numbers	Attaching Hardware Part Numbers
1035907	1086112-NLA
(2) Screw, PPH (1/4-20x2-1/2") (2) Spacer (1/4x11/32x1-3/4") (2) Locknut (1/4-20)	(2) Screw, PPH (1/4-20x2-1/2") (2) Spacer (1/4x11/32x1-3/4") (2) Locknut (1/4-20)

Footrest and Legrest Accessories

AT5543 Heel Loops and Hardware

Heel Loops and Hardware AT5543 Adjustable Angle	Heel Loops with Ankle Strap and Hardware AT5543 Adjustable Angle
1160723	1160724
(2) Heel Loops (4) Screw, PPH (1/4-20 x 2-1/4") (4) Spacer (4) Locknut (1/4-20)	(2) Heel Loops w/ Ankle Straps (4) Screw, PPH (1/4-20 x 2-1/4") (4) Spacer (4) Locknut (1/4-20)

Toe Loops (1345) and Attaching Hardware

Footplate Type	Toe Loop w/o Hardware Single #	Toe Loop w/o Hardware Pair #	Attaching Hardware Part Number	Attaching Hardware Consists of
1651, 1650, 1450, 1350, Footplates	04572X025	1037087	1037089	(4) Screw, Phillips Pan Head (10-32 x 7/8") (4) Washer (1/4 x 3/4 x 3/32") (4) Locknut (10-32)
AT5543 Footplates	04572X025	1037087	1104749	(4) Screw, Phillips Pan Head (10-32 x 1/2") (4) Washer (1/4 x 3/4 x 3/32")