




Heel Loops and Ankle Straps



1651 Composite

	Heel Loop Single #	Heel Loop Pair #	Heel Loop w/ Ankle Strap Single # - Right	Heel Loop w/ Ankle Strap Single # - Left	Heel Loop w/ Ankle Strap Pair #	Attaching Hardware Pair #
Small	1025234-NLA	1035905	1035908-NLA	1035909-NLA	1035906	1035931
Medium	1017621	1035917-NLA	1053313-NLA	1053314-NLA	1053329	1035931

Heel Loops and Ankle Straps

Adjustable Angle 5542 and AT5543

Footplate Width	Heel Loop Single #	Heel Loop Pair #	Heel Loop w/ Ankle Strap Single # - Right	Heel Loop w/ Ankle Strap Single # - Left	Heel Loop w/ Ankle Strap Pair #	Attaching Hardware Pair #
3-1/4" Wide Footplate	1018269	1036099-NLA	1036097-NLA	1036097-NLA	1036102	1035931
4-1/4" Wide Footplate	1018270	1036100	1027182	1027182	1036103-NLA	1035931
5 & 6-1/4" Wide Footplate	1030958-NLA	1036101	1036098	1036098	1036104-NLA	1035931
7-1/4" Wide Footplate	1069037	N/A	1069049-NLA	1069049-NLA	N/A	1035931

Heel and Toe Loop Attaching Hardware Kits

1035931 Kit, Heel Loop Attaching Hardware (Pair)	1037089 Hardware, Toe Loop (Pair)
Kit Includes:	Kit Includes:
(2) Screw, Phillips, Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Hex Nut w/ Nylon Cap (1/4-20) (2) Flat Washer (1/4 x 5/8 x 1/16")	(4) Screw, Phillips, Pan Head (10-32 x 7/8") (4) Washer (1/4 x 3/4 x 3/32") (4) Locknut (10-32)