



AT5543 Adjustable Angle Footplates

60 deg., 60 deg. MFX, 70 deg., 70 deg. MFX, PW93, PW93E, and 90 deg. Swingaway Supports



AT5543 Adjustable Angle Footplates

60 deg., 60 deg. MFX, 70 deg., 70 deg. MFX, PW93, PW93E, and 90 deg. Swingaway Supports

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Pivot/Slide Tube and Footplate Assemblies	1	1
2	B,F	1113572	Footplate, Adjustable Angle, Medium, Black (4-1/4" W x 8" L)	1	1
2	C,F	1113573	Footplate, Adjustable Angle, Large, Black (5-1/4" W x 8" L)	1	1
2	D,F	1113574	Footplate, Adjustable Angle, Extra Large, Black (6-1/4" W x 8" L)	1	1
3	E	1058517	Clamp, Footplate	1	1
4		1057855	Tube, Pivot and Slide, Straight, Steel, Black (9-7/8" L)	1	1
4		1057856	Tube, Pivot/Slide, Straight, Steel, Black (11-7/8" L)	1	1
5		1048816	Screw, Nylon Adjustment (1/2-13 x 7/8")	1	1
6		1052102	Hinge, Footplate - Right	1	1
6		1052103	Hinge, Footplate - Left	1	1
7	E	1029545	Screw, Socket Flat Head (1/4-20 x 3/4")	1	2
8		7930321-NLA	NLA Package, Screw, Socket Button Head (1/4-20 x 1-7/8")	25	1
9		1026900	Washer (1/4 x 23/32 x 1/16")	1	1
10		7930313	Package, Washer, Nylon (5/16 x 7/8 x 1/8")	25	1
11		1036435	Washer, Nylon (1/4 x 7/8 x 1/32")	1	1
12		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	1
13		1026099	Spacer, Coved (1/4 x 11/16 x 1/2")	1	1
14		1026907	Washer (1/4 x 1/2 x 1/16")	1	1
15		7930333	Locknut (1/4-20)	25	1
16	E	1082686	Nut, Barrel, Hex Head (1/4-20 x 11/16")	1	2

NOTE: A - Includes Item 2 thru 16.
 B - For 12-13" wide chairs.
 C - For 14-15" wide chairs.
 D - For 16" and wider chairs.
 E - Items 3, 7 and 16 must be ordered together for chairs manufactured before 8/10/99 due to a hardware change.
 F - This footplate was revised 9/2/2008 with rounder corners and accessory holes changed from 8 to 4. Recommended to replace footplates in pairs to match sides.

#1 Pivot/Slide Tube and Footplate Assemblies

Rigging Type	90, 70 MFX w/o Pivot/Slide Tube	60, 60 MFX, 70 PW93, PW93E (9-7/8" Tube)	60, 60 MFX, 70 PW93, PW93E (11-7/8" Tube)
12-13" W - Right (4-1/4" Footplate)	1060436-NLA	1060442	N/A
12-13" W - Left (4-1/4" Footplate)	1060437	1060443	N/A
14-15" W - Right (5-1/4" Footplate)	1060438	1060444	1062758
14-15" W - Left (5-1/4" Footplate)	1060439	1060445	1062759

AT5543 Adjustable Angle Footplates

60 deg., 60 deg. MFX, 70 deg., 70 deg. MFX, PW93, PW93E, and 90 deg.
Swingaway Supports

#1 Pivot/Slide Tube and Footplate Assemblies

Rigging Type	90, 70 MFX w/o Pivot/Slide Tube	60, 60 MFX, 70 PW93, PW93E (9-7/8" Tube)	60, 60 MFX, 70 PW93, PW93E (11-7/8" Tube)
16" and Wider - Right (6-1/4" Footplate)	1062570	1062756	1062760
16" and Wider - Left (6-1/4" Footplate)	1062571	1062757	1062761