

Footrest and Legrest Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	B	1027644-NLA	NLA - Strap, Calf w/ Hook and Loop (26" L)	1	1
	C	1032243	Strap, Calf w/ Hook and Loop (36" L)	1	1
	A	8881038231	Kit, 3" Extension w/ Hardware (Pair)	2	1
	D	See Chart	Footplate Type	1	1
	D	See Chart	AT5543 Heel Loops and Hardware	1	1
	D	See Chart	Toe Loops and Attaching Hardware	1	1

NOTE: A - Specify frame color when ordering.
 B - 10-14" Wide
 C - 14-16" Wide
 D - Not Shown

Footplate Type

Composite (1651) Footplate	Heel Loop Single #	Heel Loop Pair #	Heel Loop w/ Ankle Strap Right - Single #	Heel Loop w/ Ankle Strap Left - Single #	Heel Loop w/ Ankle Strap Pair #	Attaching Hardware Pair #
Small 13-14" Wide	1028412-NLA	1035912-NLA	1035908-NLA	1035909-NLA	1035906	1035931
Medium 15-16" Wide	1017621	1035917-NLA	1035919	1035920	1035918	1035931

AT5543 Heel Loops and Hardware

Footplate Type	Part Number	Consists of:
AT5543 Footplate	1160723 Kit, Heel Loop w/ Hardware, Adjustable Angle Footplate (Pair)	(2) Heel Loops (4) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (4) Spacer (4) Locknut (1/4-20)
AT5543 Footplate	1160724 Kit, Heel Loop w/ Ankle Strap and Hardware, Adjustable Angle Footplate (Pair)	(2) Heel Loops w/ Ankle Straps (4) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (4) Spacer (4) Locknut (1/4-20)
Attaching Hardware Only	1035931 Kit, Heel Loop Attaching Hardware (Pair)	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Hex Nut w/ Nylon Cap (1/4-20) (2) Flat Washer (1/4 x 5/8 x 1/16")

Footrest and Legrest Accessories

Toe Loops and Attaching Hardware

Toe Loops (1345)	Toe Loop Single #	Toe Loop Pair #
All Sizes	04572X025	1037087
Footplate Attaching Hardware	Part Number	Consists of:
1651 Footplate	1037089 Hardware, Toe Loop (Pair)	(4) Screw, Phillips Pan Head (10-32 x 7/8") (4) Washer (1/4 x 3/4 x 3/32") (4) Locknut (10-32)
AT5043 and AT5543 Footplates	1104749 Hardware, Toe Loops and Footplate (Pair)	(4) Screw, Phillips Pan Head (10-32 x 1/2") (4) Washer (1/4 x 3/4 x 3/32")