


**Pivot Slide Tube with Footplates**  
**For 93E, 93M, 93M2, H93E, H93M, H93M2 Assemblies**  
**(16" to 24" Chairs)**



# Pivot Slide Tube with Footplates

## For 93E, 93M, 93M2, H93E, H93M, H93M2 Assemblies

### (16" to 24" Chairs)

| Item Number                                                                                                                                                                                                | Note | Part Number | Description                                                              | Quantity Per |          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|-------------|--------------------------------------------------------------------------|--------------|----------|
|                                                                                                                                                                                                            |      |             |                                                                          | Package      | Assembly |
| 1                                                                                                                                                                                                          | A    | 7500084-NLA | NLA - Pivot/Slide Tube w/ Plug Button, Chrome                            | 1            | 1        |
| 2                                                                                                                                                                                                          |      | 1027097     | Plug Button, Black (3/4")                                                | 1            | 2        |
| 3                                                                                                                                                                                                          |      | 1044213-NLA | NLA - Footplate with Rubber Edge Assembly (Pair)                         | 2            | 1        |
| 3                                                                                                                                                                                                          |      | 1043659-NLA | NLA - Footplate, Aluminum, Deluxe, Black - Right (8" W x 5-3/4" L)       | 1            | 1        |
| 3                                                                                                                                                                                                          |      | 1043660-NLA | NLA - Footplate, Aluminum, Deluxe, Black - Left (8" W x 5-3/4" L)        | 1            | 1        |
| 4                                                                                                                                                                                                          | A,B  | 7930354-NLA | NLA - Package, Phillips Head Screw (#8-32 x 3/8")                        | 25           | 1        |
| 5                                                                                                                                                                                                          |      | 1036120-NLA | NLA - Spring                                                             | 1            | 2        |
|                                                                                                                                                                                                            |      | 1028905-NLA | NLA - Pivot/Slide Tube w/ Rubber Edge Footplate Assembly, Chrome - Right | 1            | 1        |
|                                                                                                                                                                                                            |      | 1028906-NLA | NLA - Pivot/Slide Tube w/ Rubber Edge Footplate Assembly, Chrome - Left  | 1            | 1        |
| NOTE: A - Front Riggings built before 12/93 had footplates with a flat black finish. After 11/93 the finish is "Wrinkle Black." For Front Riggings before 12/93, order in pairs to match.<br>B - Not Shown |      |             |                                                                          |              |          |