






MK5 SPJ (After 6/26/03) & SPJ+ (Introduced 8/16/05) Joystick and Mounting Hardware

The SPJ+ Joystick plus Cable Kit replaces the SPJ Joystick after 10/31/07.



ASBA Seat, Fixed Height Joystick Offset Tube Clamp Assembly

MK5 SPJ (After 6/26/03) & SPJ+ (Introduced 8/16/05) Joystick and Mounting Hardware

The SPJ+ Joystick plus Cable Kit replaces the SPJ Joystick after 10/31/07.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,G	1153650	Kit, SPJ+ Joystick w/ Cable	1	1
2			Joystick, SPJ+	1	1
3			Knob, Joystick Rubber, Black	1	1
4			Skirt, Joystick Rubber, Black	1	1
5			NLA - Overlay, SPJ+ Joystick	1	1
6			Cable, Joystick SPJ+	1	1
7		1116404	Cable, Joystick Extension	1	1
2		1127291	Joystick, SPJ+	1	1
	H	See Note-NLA	See Footnote	1	1
3		1040217	Knob, Joystick Rubber, Black	1	1
4		1040218	Skirt, Joystick Rubber, Black	1	1
5		1165094	Kit, Overlay, SPJ+ Joystick (MK5/MK6)	1	1
			NLA - Overlay, SPJ+ Joystick	1	1
			Instruction Sheet, MK5/MK6 Overlay	1	1
6		1127293	Cable, Joystick SPJ+	1	1
7		1116404	Cable, Joystick Extension	1	1
8	E	1161733	Kit, 2" Offset Joystick Tube w/ Hardware	1	1
9			Plug Button, Black (7/8")	1	1
10			Tie Wrap (11-1/2" L)	1	2
11			Screw, Phillips Pan Head (#10-32 x 1-1/8")	1	2
11			Screw, Phillips Pan Head (M5 x .8 x 30mm)	1	2
12			Spacer, Joystick Tube, Black	1	1
			Tube, Joystick Offset (2")	1	1
8	B	1093069	Kit, 1/2" Offset Joystick Tube with Hardware	1	1
9		1027095	Plug Button, Black (7/8")	1	1
10		1095502-NLA	NLA - Package, Tie Wrap (11-1/2" L)	10	1
11		1025714	Screw, Phillips Pan Head (#10-32 x 1-1/8")	1	2
12		1093063	Spacer, Joystick Tube, Black	1	1
13	C	1115000	Kit, Joystick Tube Mounting Hardware - Captain Seat	1	1
14		1105049	Clamp Assembly, Fixed Height, Captain/ Van Seat	1	1
15		1037610	Screw, Hex Head (1/4-20 x 1-5/8")	1	3
16		1087219	Spacer (1/4 x 3/8 x 1")	1	3
17		1025195	NLA - Locknut (1/4-20)	1	3
18		1139422	Clamp, Height Adjustable Joystick Tube (ARM250)	1	1
19		1093071	Clamp, Fixed Height Joystick Tube	1	1
	D	1006270	Extension, Joystick, "T" Style (1560) (4" W x 5-1/2" L)	1	1
	D	1006271	Extension, Joystick, Straight Style (1561) (1/2" W x 5" L)	1	1

MK5 SPJ (After 6/26/03) & SPJ+ (Introduced 8/16/05) Joystick and Mounting Hardware

The SPJ+ Joystick plus Cable Kit replaces the SPJ Joystick after 10/31/07.

NOTE: A - Includes items 2 thru 6.
B - Used on A.S.B.A. seats. Includes Joystick Tube plus items 9 thru 12.
C - Includes items 14 thru 17.
D - Not Shown.
E - Used on Van Seats. Includes Joystick Tube and items 9 thru 13. Kit contains Standard (#10-32) and Metric (M5) screws. Customer selects the screw that applies to their application.
G - Before replacing your joystick, verify with Invacare Technical Service Phone Support first to see if only a Joystick Overlay Kit is needed
H - MK5 SPJ Joystick part number is no longer available as a service part. Use item #1 1153650 as the replacement.