


















# Electronics - Arro BA2180



R2009\_082\_002A

## Electronics - Arro BA2180

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		MTRM25031	Actuator, Head (Red) Linak # LA31-U086	1	1
2		MTRM25021	Actuator, Foot (Blue) Linak # LA31-U087	1	1
3	B	MTRM25471-NL A	NLA - Actuator, Hilo, Foot, Short Cord (White) Linak # LA31-U010	1	1
4		MTRM25472	Actuator, Hilo, Head, Long Cord (Yellow) Linak # LA31-U011	1	1
5	B,D	JNCTNDC1B	Control Box (111) Linak # CB09-u011	1	1
5	C,D	JNCTNDC1D	Box, Control (111) Linak # CB09-U044-XX	1	1
5	D	JNCTNDC1C-NL A	NLA - Box, Control (151/171) Linak # CB09-U037	1	1
6		JNCTNDC35	Box, DJB (151/171) Linak # DJBX-001-XX	1	1
7		1190227	Pendant, Six Function (w/ Loop Adhesive Back Strip and DB9 Connector) Linak # PNDNT-6FNC-1-2	1	1
8		1190232	Pendant, Six Function (w/ Strap and DB9 Connector) Linak # PNDNT-6FNC-1-2	1	1
9		PNDNTYCABLE	Cable, Y, Pendant	1	1
10		1190225	Cable, Conversion	1	1
11		PNDNT6FNCIV3	Pendant, Six Function w/ RJ48 Connector) Linak # PNDNT-6FNC-1-2	1	1
12	E	JNCTNDC34-NL A	NLA - Cable, T, Pendant, Fine Thread Nut (151/171)	1	1
12	B	JNCTNDC47A-N LA	NLA - Cable, T, Pendant, Coarse Thread Nut (151/171)	--	--
13		JNCTNDC202	Pad, Membrane (151)	1	1
14		JNCTNDC20A2	Pad, Membrane w/ CPR (171)	1	1
15	F	JNCTNDC20D-N LA	NLA - Pad, Membrane w/o CPR (171 w/o CPR)	1	1
16		JNCTNBOXKIT3 -NLA	NLA - Box, Lockout (111)	1	1
17		1152564	Battery, On Board Backup w/ Hardware	1	1
18		JNCTNDCFUUSE K2-NLA	NLA - Kit, Fuse, Arro, 2.5 Amps (151/171)	1	1
18		JNCTNDCFUUSE K1	Kit, Fuse (1.6 Amps)	1	1
19		JNCTNDC31-NL A	NLA - Cable, Battery	1	1
20		JNCTNDCCORD KI	Cable, Power	1	1
21		JNCTNDC32	Cord, Membrane (Footboard)	1	1
22		JNCTNDC33-NL A	NLA - Cable, Membrane (DJB)	1	1

NOTE: B - Used on Beds made Before 8/26/04. Use MTRM25472  
C - Used on Beds made After 8/25/04.  
D - Control Box type determined by the item number on the box. In the item number CB09-Uxxx-(xx), the last two digits are a revision level and have nothing to do with part identification.  
E - Order JNCTNDC47A (Coarse Thread) to match Pendant thread.  
F - Order JNCTNDC20A2 (Pad, Membrane)