



# GB Connector Covers



## GB Connector Covers

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1141670-NLA	NLA - GB Connector Cover w/o Wheel Lock Assembly - Right	1	1
2			Bracket, Cover, Connector, Black - Right	1	1
3			Bracket, Cover, Connector, Top, Black - Right	1	1
4			NLA - Screw, Phillips Pan Head w/ Patch (#10-32 x 3/8")	1	4
5			NLA - Washer, Plain (3/16 x 1/2 x 1/16")	1	2
			Sheet, Instruction	1	1
1	A	1141671-NLA	NLA - GB Connector Cover w/o Wheel Lock Assembly - Left	1	1
2			Bracket, Cover, Connector, Black - Left	1	1
3			Bracket, Cover, Connector, Top, Black - Left	1	1
4			NLA - Screw, Phillips Pan Head w/ Patch (#10-32 x 3/8")	1	4
5			NLA - Washer, Plain (3/16 x 1/2 x 1/16")	1	2
			Sheet, Instruction	1	1
6	B	1141672	GB Connector Cover w/ Wheel Lock Assembly - Right	1	1
2			Bracket, Cover, Connector, Black - Right	1	1
4			NLA - Screw, Phillips Pan Head w/ Patch (#10-32 x 3/8")	1	2
5			NLA - Washer, Plain (3/16 x 1/2 x 1/16")	1	2
			Sheet, Instruction	1	1
6	B	1141673	GB Connector Cover w/ Wheel Lock Assembly - Left	1	1
2			Bracket, Cover, Connector, Black - Left	1	1
4			NLA - Screw, Phillips Pan Head w/ Patch (#10-32 x 3/8")	1	2
5			NLA - Washer, Plain (3/16 x 1/2 x 1/16")	1	2
			Sheet, Instruction	1	1

NOTE: A - Includes items 2-5 and instruction sheet.  
 B - Includes items 2, 4, 5 and instruction sheet.