


Pivot Slide Tube with Footplate For AL4A, AHL4A Assemblies (16" to 24" Chairs)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1016346	Cap, End, Tube, Black (7/8")	1	1
2	A,C	44200X027	Footplate, Aluminum, Medium, Black (7-3/4" W x 6" L)	1	1
3		7930354-NLA	NLA - Package, Phillips Head Screw (#8-32 x 3/8")	1	1
4		1006981	Spring, Tension	25	2
5		1041526	Pivot/Slide Tube w/ Plug Button Assembly, Chrome	1	1
6		1027097	Plug Button, Black (3/4")	1	1
7		1018133	Plug Button, Black (5/8")	1	1
8		1070168	Plug, Push, Footplate	1	1
	A,B	1023907	Footplate w/Pivot Slide Tube - RH (18")	1	1
	A,B	1023908	Footplate w/Pivot Slide Tube - LH (18")	1	1

NOTE: A - Aluminum footplates on chairs built before 2/96 had a "Bright Aluminum" finish. Footplates after 1/96 now have a black finish. For Front Riggings before 2/96, order in pairs to match.

B - Not Shown

C - Right or Left