



Side Rail Parts - BA0780


Assist Rails




Side Rails - Non Tuck Including Hardware



Rail Arms



Pull Knob Assembly for All Rails



Side Rail Parts - BA0780

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		RAILSASSTAEP R	Assist Rails (Pair), "U" Bracket Style	1	1
2		RAILSASSTAEL T-NLA	NLA - Assist Rail - Left Side	1	1
3		RAILSASSTAER T-NLA	NLA - Assist Rail - Right Side	1	1
4		KITASSISTRAIL	Kit, Assist Rail Mounting Hardware	1	1
5		RLBRACKASSIS T-NLA	NLA - U Bracket, Three Position - Left Only	1	1
6		RRBRACKASSI ST-NLA	NLA - U Bracket, Three Position - Right Only	1	1
7		RLBRACKASTR UP-NLA	NLA - U Bracket, No Guard - Left Hand Only	1	1
8		RRBRACKASTR UP-NLA	NLA - U Bracket, No Guard - Right Hand Only	1	1
9		IHARRLFUL	NLA - Full Rails (Set)	1	1
10		RAILS0080FLNT -NLA	NLA - Full Rail - Left Hand	1	1
11		RAILS0080FRN T-NLA	NLA - Full Rail - Right Hand	1	1
12		IHARRLPAR	NLA - Partial Rails (Set)	1	1
13		RAILS0080PLNT -NLA	NLA - Partial Rail - Left Hand	1	1
14		RAILS0080PRN T-NLA	NLA - Partial Rail - Right Hand	1	1
15		LRA001RA-NLA	NLA - Long Rail Arm - Right	1	1
16		SRA001RA-NLA	NLA - Short Rail Arm - Right	1	1
17		SRA001LA-NLA	NLA - Short Rail Arm - Left	1	1
18		LRA001LA-NLA	NLA - Long Rail Arm - Left	1	1
19		RAILBRACKNT- NLA	NLA - Rail Bracket	1	1
20		COL34R-NLA	NLA - Locking Set Collar	1	1
21		1189259	Kit, Knob - Latch	1	1
			Knob, LATCH-0002 (3/8-16 x 1-3/8")	1	1
			Spring, Compression (3/4 x 1 x 3/64")	1	1
			Screw, Shoulder Slotted Head (1/2 x 1/2", 3/8-16)	1	1
22		1186856	Spring, Compression (3/4 x 1 x 3/64")	1	2
23		1165068	Package, End Plug	1	1
24		1189260	Kit, Assist Rail Mounting	1	1
			Bolt, Shoulder (3/8-16 x 3/4")	1	1
			Washer, Belleville (33/64 x 1-1/64 x 3/32")	1	1
			Washer, Flat Edge (17-32 x 1-1/4 x 1/8")	1	1
			Washer, Nylon (33/64 x 2-1/2 x 3/64")	1	1
			Locknut (3/8-16)	1	1
			Washer, Plain (17/32 x 1-1/4 x 1/8")	1	1
25		1189261	Kit, Rail Mounting Hardware	1	1
			Pin, Rail (5/16-18)	1	4
			Washer, Nylon (33/64 x 1-29/64 x 3/64")	1	4
			Washer, Retaining (11/32 x 1-3/32 x 3/64")	1	4

Side Rail Parts - BA0780

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
			NLA - Nut, Hex Head, Flange (5/16-18)	1	4
			Cap, Bearing	1	4