

Spoke Rear Wheels (w/o Handrims) (Before 9/2/11)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	1030936	Plug, Valve Stem	1	1
	A	See Chart	Wheels without Handrims - Spoke	1	1

NOTE: A - Not Shown

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
20SP	20 x 1-3/8"	Pneumatic	1115585-NLA	1017519	1017520	40117X021	1015756
20SP	20 x 1-3/8"	Pneumatic w/ Flat Free	1115589-NLA See Note(s) A	1017519	1023322	40117X021	1015756
22SP	22 x 1-3/8"	Pneumatic	1115656-NLA	1018544	1018545	40117X021	1015756
22SP	22 x 1-3/8"	Pneumatic w/ Flat Free	1115658-NLA See Note(s) A	1018544	1023323	40117X021	1015756
24SP	24 x 1-3/8"	Pneumatic	1115660	40176X060	40176X061	40117X021	1015756
24SP	24 x 1-3/8"	Pneumatic w/ Flat Free	1116707 See Note(s) A	40176X060	1010906	40117X021	1015756
24SP	24 x 1"	Primo	1115662-NLA	1026138	1032240	40117X021	1015756
24SP	24 x 1"	KIK	1114112	1062852	N/A	40117X021	1015756
24SP	24 x 1"	Schwalbe	1154227	1139813	1105941	40117X021	N/A
24HPS	24 x 1-3/8"	Pneumatic	1116709-NLA See Note(s) E	40176X060	40176X061	40117X021	1015756
24HPS	24 x 1-3/8"	Pneumatic (Pair)	1116710	40176X060	40176X061	40117X021	1015756
24HPS	24 x 1-3/8"	Pneumatic w/ Flat Free	1116715 See Note(s) A	40176X060	1010906	40117X021	1015756
24HPS	24 x 1"	Primo	1116711	1026138	1032240	40117X021	1015756
24HPS	24 x 1"	KIK	1114138	1062852	N/A	40117X021	1015756
24HPS	24 x 1"	Schwalbe	1157361-NLA	1139813	1105941	40117X021	N/A
SPIN24 Black Spokes	24 x 1"	Primo	1112338-NLA See Note(s) C	1026138	1105941	1104350	N/A
SPIN24 Yellow, Red, Blue or White Spokes	24 x 1"	Primo	8881159332 See Note(s) D	1026138	1105941	1104350	N/A
SPIN24 Black Spokes	24 x 1"	Schwalbe	1154229	1139813	1105941	1104350	N/A
SPIN24 Yellow, Red, Blue or White Spokes	24 x 1"	Schwalbe	8881159332 See Note(s) D	1139813	1105941	1104350	N/A
SPLX24 Black, Yellow, Red, Blue or White Spokes	24 x 1"	Primo	8881159332 See Note(s) D	1026138	1105941	1104350	N/A

Spoke Rear Wheels (w/o Handrims) (Before 9/2/11)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
SPLX24 Black, Yellow, Red, Blue or White Spokes	24 x 1"	Schwalbe	8881159332 See Note(s) D	1139813	1105941	1104350	N/A
FN1624	24 x 1"	Primo	8881159332 See Note(s) D	1026138	1105941	40117X021	N/A
FN1624	24 x 1"	Pneumatic	8881159332 See Note(s) D	40176X060	40176X061	40117X021	N/A
FN1624	24 x 1"	KIK	8881159332 See Note(s) D	1062852	N/A	40117X021	N/A
FN1624	24 x 1"	Schwalbe	8881159332 See Note(s) D	1139813	1105941	40117X021	N/A
559TS	25 x 1"	Primo	1116717-NLA	1044947	1029471	40117X021	1015756
599TS	25 x 1"	Primo Treaded Pneumatic	1130592-NLA	1130527	1104326	40117X021	1015756
SPIN25 Black Spokes	25 x 1"	Primo	1112348-NLA See Note(s) C	1044947	1104326	1104350	N/A
SPIN25 Yellow, Red, Blue or White Spokes	25 x 1"	Primo	8881159332 See Note(s) D	1044947	1104326	1104350	N/A
SPIN25 Black Spokes	25 x 1"	Primo Treaded Pneumatic	1130593-NLA See Note(s) C	1130527	1104326	1104350	N/A
SPIN25 Yellow, Red, Blue or White Spokes	25 x 1"	Primo Treaded Pneumatic	8881159332 See Note(s) D	1130527	1104326	1104350	N/A
SPIN25 Black Spokes	25 X 1"	Schwalbe	1154232-NLA See Note(s) B	1139814	1104326	1104350	N/A
SPIN25 Yellow, Red, Blue or White Spokes	25 x 1"	Schwalbe	8881159332 See Note(s) D	1139814	1104326	1104350	N/A
SPLX25 Black, Yellow, Red, Blue or White Spokes	25 x 1"	Primo	8881159332 See Note(s) D	1044947	1104326	1104350	N/A
SPLX25 Black, Yellow, Red, Blue or White Spokes	25 x 1"	Primo Treaded Pneumatic	8881159332 See Note(s) D	1130527	1104326	1104350	N/A
SPLX25 Black, Yellow, Red, Blue or White Spokes	25 x 1"	Schwalbe	8881159332 See Note(s) D	1139814	1104326	1104350	N/A
FN1625	25 x 1"	Primo	8881159332 See Note(s) D	1044947	1104326	40117X021	N/A

Spoke Rear Wheels (w/o Handrims) (Before 9/2/11)

Wheel Style	Wheel Size	Tire Type	Wheel Assembly	Tire	Tube/Insert	Bearing	Hub Spacer
FN1625	25 x 1"	Primo Treaded Pneumatic	8881159332 See Note(s) D	1130527	1104326	40117X021	N/A
FN1625	25 x 1"	KIK	8881159332 See Note(s) D	1144435	N/A	40117X021	N/A
FN1625	25 x 1"	Schwalbe	8881159332 See Note(s) D	1139814	1104326	40117X021	N/A
26HPS	26 x 1-3/8"	Pneumatic	1115652-NLA	1026142	1026140	40117X021	1015756
26HPS	26 x 1-3/8"	Pneumatic	1116725-NLA	N/A	N/A	N/A	N/A
26HPS	26 x 1"	Primo	1115654	1029470	1029471	40117X021	1015756
26HPS	26 x 1"	Schwalbe	1154389-NLA	1139815	1104326	40117X021	N/A
SPIN26 Black Spokes	26 x 1"	Primo	1112343-NLA See Note(s) C	1029470	1104326	1104350	N/A
SPIN26 Yellow, Red, Blue or White Spokes	26 x 1"	Primo	8881159332 See Note(s) D	1029470	1104326	1104350	N/A
SPIN26 Black Spokes	26 x 1"	Schwalbe	1154235-NLA	1139815	1104326	1104350	N/A
SPIN26 Yellow, Red, Blue or White Spokes	26 x 1"	Schwalbe	8881159332 See Note(s) D	1139815	1104326	1104350	N/A
FN1626	26 x 1"	Primo	8881159332 See Note(s) D	1029470	1104326	40117X021	N/A
FN1626	26 x 1"	Pneumatic	8881159332 See Note(s) D	1026142	1026140	40117X021	N/A
FN1626	26 x 1"	Schwalbe	8881159332 See Note(s) D	1139815	1104326	40117X021	N/A
26SP	26 x 1-3/8"	Pneumatic	1116719-NLA	1026142	1026140	40117X021	1015756
26SP	26 x 1"	Primo	1116723-NLA	1029470	1029471	40117X021	1015756
26SP	26 x 1"	Schwalbe	1154390-NLA	1139815	1104326	40117X021	N/A

NOTE: A - Stem Hole Plugs (Part Number 1030936) are used with Flat Free and Urethane Tires.
B - Tire is low profile and dimension stated includes overall diameter.
C - As of 11/17/08 Spinergy Wheels changed from black spoke, black hub to a black spoke, silver hub appearance. For chairs built Before 11/17/08, it is recommended when one wheel is requested that the wheels be replaced in pairs.
D - Use the Parts Wheel Configurator when ordering a complete wheel for this model. To omit handrim in the Parts Wheel Configurator, order as a special
E - Wheels must be replaced in pairs. See Part number 1116710 for replacement.