

# Handrims and Mounting Hardware

For Natural Fit Handrims, please refer to that particular category.

## Handrims

| Wheel Style | Anodized Aluminum            | Plastic Coated<br>(Black) | Titanium | Attachment<br>Points |
|-------------|------------------------------|---------------------------|----------|----------------------|
| 20COM       | 1033374                      | 1115400                   | N/A      | 7                    |
| 20SP        | 1043442-NLA                  | 1043688                   | N/A      | 4                    |
| 20ITESPK    | 1043442-NLA                  | 1043688                   | N/A      | 4                    |
| 22COM       | 1033373                      | 1115401                   | N/A      | 7                    |
| 22SP        | 1043443                      | 1043689-NLA               | N/A      | 4                    |
| 22ITESPK    | 1043443                      | 1043689-NLA               | N/A      | 4                    |
| 24COM7      | 1033372                      | 1115397                   | N/A      | 7                    |
| 24SP        | 1026136                      | 1035612                   | 1060873  | 6                    |
| 24ITESPK    | 1026136                      | 1035612                   | 1060873  | 6                    |
| SPIN24      | 1114300                      | N/A                       | 1060873  | 6                    |
| SPLX24      | 1114300                      | N/A                       | 1060873  | 6                    |
| FN1624      | 1114300                      | N/A                       | 1060873  | 6                    |
| 24HPS       | 1026136                      | 1035612                   | 1060873  | 6                    |
| 559TS       | 1029493                      | 1043249                   | N/A      | 6                    |
| SPIN25      | 1114301-NLA                  | N/A                       | 1060874  | 6                    |
| SPLX25      | 1114301-NLA                  | N/A                       | 1060874  | 6                    |
| FN1625      | 1114301-NLA                  | N/A                       | 1060874  | 6                    |
| 26SP        | 1026137                      | 1035613                   | 1060875  | 6                    |
| 26HPS       | 1026137                      | 1035613                   | 1060875  | 6                    |
| SPIN26      | 1114302-NLA<br>See Note(s) A | N/A                       | 1060875  | 6                    |
| FN1626      | 1114302-NLA<br>See Note(s) A | N/A                       | 1060875  | 6                    |

NOTE: A - Part No Longer Available individually. See Part # 1112344 for replacement.

# Handrims and Mounting Hardware

For Natural Fit Handrims, please refer to that particular category.

## Handrim Hardware

| Wheel Style | Anodized Aluminum | Plastic Coated<br>(Black) | Titanium |
|-------------|-------------------|---------------------------|----------|
| 20COM       | 1025153           | 1025158                   | N/A      |
| 20SP        | 1043708           | 1043708                   | N/A      |
| 20ITESPK    | 1043708           | 1043708                   | N/A      |
| 22COM       | 1026384           | 1026383                   | N/A      |
| 22SP        | 1043708           | 1043708                   | N/A      |
| 22ITESPK    | 1043708           | 1043708                   | N/A      |
| 24COM7      | 1026384           | 1026383                   | N/A      |
| 24SP        | 1026910           | 1026910                   | 1026910  |
| 24ITESPK    | 1026910           | 1026910                   | 1026910  |
| SPIN24      | 1026910           | N/A                       | 1026910  |
| SPLX24      | 1026910           | N/A                       | 1026910  |
| FN1624      | 1026910           | N/A                       | 1026910  |
| 24HPS       | 1026910           | 1026910                   | 1026910  |
| 559TS       | 1026910           | 1026910                   | N/A      |
| SPIN25      | 1026910           | N/A                       | 1026910  |
| SPLX25      | 1026910           | N/A                       | 1026910  |
| FN1625      | 1026910           | N/A                       | 1026910  |
| 26SP        | 1026910           | 1026910                   | 1026910  |
| 26HPS       | 1026910           | 1026910                   | 1026910  |
| SPIN26      | 1026910           | N/A                       | 1026910  |
| FN1626      | 1026910           | N/A                       | 1026910  |