


Miscellaneous Accessories

Pelvic Belt is not available as a service part. Transport Ready chairs must be returned to Invacare for service with regard to specific TRRO parts.



Auto Style Push Button Buckle Positioning Strap



Hook and Loop Style Buckle Positioning Strap

Miscellaneous Accessories

Pelvic Belt is not available as a service part. Transport Ready chairs must be returned to Invacare for service with regard to specific TRRO parts.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Seat/Chest Positioning Straps	1	1
	A,B	1052425	Kit, Chest Positioning Strap Mounting Hardware	2	1
	A	1008495	Screw, Socket Button Head (1/4-20 x 1-1/2")	1	2
	A	1025736	Locknut (1/4-20)	1	2
		See Chart	Miscellaneous Accessories	--	--

NOTE: A - Not Shown
B - Includes 2 ea. Screws and 2 ea. Locknuts

#1 Seat/Chest Positioning Straps

Chair Width	Push Button (Auto) Style Buckle	Hook and Loop Style Fastener
14-19" (45" Long)	1004310	1034744
20" (61" Long)	1034625	1034745

Pelvic Belt is not available as a service part. Transport Ready Chairs requiring any parts related to the Transport Ready Option must be returned to Invacare for service.

Miscellaneous Accessories

Part Number	Description	Package Quantity	Assembly Quantity
1130574 See Note(s) A	Anti-Tipper	1	2
1109322-NLA	Back Pack, 14" Wide Chairs	1	1
1109323	Back Pack, 15-22" Wide Chairs	1	1
1027644-NLA	Calf Strap, 14-15" Wide Chairs (26" Long)	1	1
1032243	Calf Strap, 16-22" Wide Chairs (36" Long)	1	1
See Caster Swivel Locks	Caster Swivel Locks	1	1
See Fabric Clothing Protector Category	Fabric Clothing Protectors	1	See Fabric Clothing Protector Category
1032215	H-Strap, 14-18" Wide Chairs	1	1
1032216	H-Strap, 19-22" Wide Chairs	1	1
1026987 See Note(s) C	Impact Guard, 8" Long	1	2
1159758	Protector, Frame, Front, Blue	1	2

Miscellaneous Accessories

Pelvic Belt is not available as a service part. Transport Ready chairs must be returned to Invacare for service with regard to specific TRRO parts.

Miscellaneous Accessories

Part Number	Description	Package Quantity	Assembly Quantity
1159756	Protector, Frame, Back, Blue	1	1
1159384-NLA	Protector, Frame, Front, Silver	1	2
1159383-NLA	Protector, Frame, Back, Silver	1	1
1159759	Protector, Frame, Front, Red	1	2
1159757	Protector, Frame, Back, Red	1	1
1163422	Protector, Frame, Front, Black	1	2
1163421	Protector, Frame, Back, Black	1	1
See Positioning Strap Chart	Positioning Strap, Seat/Chest	1	See Positioning Strap Chart
1111934	Rear Frame Protector	1	1
See Rigid Side Guard Category	Rigid Side Guards	1	See Rigid Side Guard Category
1109324	Seat Pouch 8-1/2" W x 7" H x 4" D	1	1
1032146-NLA	Spoke Guard 20" Spoke/Composite Wheels	1	2
1032147-NLA	Spoke Guard 22" Spoke/Composite Wheels	1	2
1028442-NLA	Spoke Guard 24" Composite Wheels	1	2
1028441-NLA	Spoke Guard 24" Spoke Wheels	1	2
1028443-NLA	Spoke Guard 25"/26" Spoke Wheels	1	2
1130626 See Note(s) B	Sport Performance Detail Package	1	1
40176X010-NLA	Tire Pump	1	1

NOTE: A - For Anti-Tipper Bracket, see Axle Mounting Plate Hardware category.
 B - Includes Anodized Back Cane Tube Brackets (Right & Left), Rear Frame Protectors (2 ea.), and Impact Guards (2 ea.).
 C - For 70 deg. Tapered Footrest.