


Handrims and Handrim Hardware - Natural Fit and Natural Fit LT



6 Point Attaching Wheels (Spoke Wheels)

Handrims and Handrim Hardware - Natural Fit and Natural Fit LT

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Natural Fit Handrims - 6 Point Attaching Wheels (Spoke Wheels)	1	1
1		See Chart	#1 Natural Fit LT Handrims - 6 Point Attaching Wheels (Spoke Wheels)	1	1
2	A	1145679	Kit, Handrim Attaching Hardware - 6 Point	1	1
			Screw, Phillips Pan Head (#10-32 x 1/2")	1	6
			Locknut (#10-32)	1	6
	B	1144492	NLA - Strap, Securing	1	2
			Instruction Sheet, Natural Fit Handrims	1	1

NOTE: A - Includes Screws, Locknuts, and Securing Straps.
B - Not Shown

#1 Natural Fit Handrims - 6 Point Attaching Wheels (Spoke Wheels)

Handrims sold as 1 each. Includes all Attaching Hardware (item 2) and Instruction Sheet.

Grip Style	22" CRF6722 Wheels Only	24" Wheels	25" Wheels	26" Wheels
Standard See Note(s) A	1145545-NLA	1145546-NLA	1145547	1145548-NLA
Super Grip	1145549-NLA See Note(s) A	1145550-NLA	1145551	1145552

NOTE: A - The Natural Fit Handrims have 6 points of attachment to mount to a wheel. The 22" Spoke wheel (CRF22SP) has 4 points of attachment to mount to a wheel and will not accommodate a Natural Fit Handrim.

#1 Natural Fit LT Handrims - 6 Point Attaching Wheels (Spoke Wheels)

Handrims sold as 1 each. Includes all Attaching Hardware (item 2) and Instruction Sheet.

Grip Style	24" Wheels	25" Wheels	26" Wheels
Standard	1163646	1163647-NLA	1163648-NLA
Super Grip	1163649	1163650-NLA	1163651

NOTE: Natural Fit LT Handrims are not available in the 22" version.