


Handrims and Handrim Hardware - Natural Fit LT



7 Point Attaching Wheels (Composite Wheels)

Handrims and Handrim Hardware - Natural Fit LT

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Natural Fit LT Handrims - 7 Point Attaching Wheels (Composite Wheels)	1	1
2	B	1145680-NLA	Kit, Handrim Attaching Hardware - 7 Point	1	1
			Screw, Phillips Pan Head w/ Patch (#10-32 x 1-3/4")	1	7
			NLA - Spacer (3/16 x 5/16 x 1-1/16")	1	7
			NLA - Strap, Securing	1	2
		1144492	Instruction Sheet, Natural Fit Handrims	1	1

NOTE: B - Not Shown

#1 Natural Fit LT Handrims - 7 Point Attaching Wheels (Composite Wheels)

Handrims sold as 1 each. Includes all Attaching Hardware (item 4) and Instruction Sheet.

Grip Style	24" Wheels
Standard	1163652-NLA
Super Grip	1163653