




Canes (Walking)




**Crook with
Foam Grip Shown**



**Modified Crook with
Standard Grip shown**



**Offset Crook with
Standard Grip and Strap Shown**



**Offset with Dual Cam Grip
Shown**

Canes (Walking)

Grip Style Chart

Model	Description	Grip	Grip Kits (Includes 10 Each)	Cane Color
8901	Custom Fit Crook - Foam	1032246 See Note(s) A	1130775-NLA	Silver
8902	Custom Fit Offset - Foam	1032246 See Note(s) A	1130775-NLA	Silver
8913	Standard Crook - Foam	1032246 See Note(s) A	1130775-NLA	Black
40915-6 (Discontinued)	Crook - Standard	1048091-NLA See Note(s) A	1130776-NLA	Silver
8916	Crook - Standard	1048091-NLA See Note(s) A	1130776-NLA	Black
8917	Offset Crook- w/Strap	1048094-NLA See Note(s) A	1130777-NLA	Silver
8931	Offset - Dual Cam	1062839-NLA See Note(s) A	1130778-NLA	Silver
8932	Offset - Dual Cam	1062839-NLA See Note(s) A	1130778-NLA	Black
40918-6 (Discontinued)	Offset - Dual Cam	1062839-NLA See Note(s) A	1130778-NLA	Black
8925	Modified Crook - Standard	1048091-NLA See Note(s) A	1130776-NLA	Silver
8926	Modified Crook - Foam	1032246 See Note(s) A	1130775-NLA	Black
8962	Crook - Variety Pack			
8964	Offset - Variety Pack			

NOTE: A - Individual items no longer available, refer to Grip Kits for replacement.

Canes (Walking)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Grip Style Chart	1	1
2		40188X009	Button, Single Snap	1	1
3		1055948	Kit, Camlock	1	1
			Lever, Cam, Black	1	1
			Housing, Cam	1	1
			Instruction Sheet, Camlock	1	1
4	B,D	1065416	Tip, Rubber, Black (3/4" Tubing)	1	1
4	A,C	1086773-NLA	NLA - Tip, Rubber, Black (3/4" Tubing)	1	1
<p>NOTE: A - Item no longer available, refer to Part #1128290. B - Fits all listed except 8901 and 8902. C - To be used on Models 8901 and 8902. D - After 10/1/08, appearance of the rubber tip has changed. New tips are functionally equivalent.</p> <p>NOTE: Items 2 and 3 will fit the Canes listed above EXCEPT Custom Fit Canes 8901 and 8902.</p>					