





Lower Tiller Assembly w/ Caster



Lower Tiller Assembly w/ Caster

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1148840-NLA	NLA - Kit, Fork w/ Hardware Includes: (1) ea. Fork, (1) ea. Padlock, (1) ea. Steering Nut, (1) ea. Plain Washer, (1) ea. Tooth Washer, (1) ea. Steering Race, and (2) ea. Set Bearing.	1	1
2		1148836-NLA	NLA - Kit, Frame Connection Included: (1) ea. Frame Hook, (1) ea. Hook Pin, (1) ea. C-Ring, (1) ea. Spring, (1) ea. Lever, & (1) ea. Lever Cap.	1	1
3		1148837	Kit, Lower Steering Hardware Includes: (1) ea. Steering Adj., (1) ea. Plain Washer, (2) ea. Nylon Nut, (1) ea. Bolt, (1) ea. Joint Bolt, (1) ea. Washer, & (1) ea. Guide Complete.	1	1
4		1148859	Kit, Front Wheel Includes: (2) ea. Nylon Nuts, (1) ea. Front Wheel Assembly, (1) ea. Plain Washer, (1) ea. Bolt, & (2) ea. Spacers.	1	1
5		1144806-NLA	NLA - Boot, Bellow, Tiller	1	1