


Adductors Mounting Hardware, MPCTMT/MPCTMT-L Style

For replacement Adductors, see section titled, "Accessories-Adductors" within the Formula™CG Powered Seating Elevate CGE Tilt CGT/CGTE Recline CGR Tilt/Recline CGTR/CGTRE parts book.



Adductors Mounting Hardware, MPCTMT/MPCTMT-L Style
For replacement Adductors, see section titled, "Accessories-Adductors"
within the Formula™CG Powered Seating Elevate CGE Tilt CGT/CGTE
Recline CGR Tilt/Recline CGTR/CGTRE parts book.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1162434	Hardware, Mounting, Adductor, MPCTMT/MPCTMT-L Front Riggings	1	1
2			Plug Button, Black (1")	1	2

NOTE: A - Includes item 2.