


# Miscellaneous Accessories

Invacare does not service any parts associated with the TRRO/TRBKTS option including belts, belt attachment hardware, tie-down brackets, etc. Chairs with the Transport Ready or Transport Bracket option must be returned to Invacare for service.



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Seat/Chest Positioning Straps	1	1
	A,B	1052425	Kit, Chest Positioning Strap Mounting Hardware	2	1
	A	1008495	Screw, Socket Button Head (1/4-20 x 1-1/2")	1	2
	A	1025736	Locknut (1/4-20)	1	2
		See Chart	Miscellaneous Accessories	--	--

NOTE: A - Not Shown  
B - Includes 2 ea. Screws and 2 ea. Locknuts

## #1 Seat/Chest Positioning Straps

Chair Width	Push Button (Auto) Style Buckle	Airline Style Buckle	Hook and Loop Style Fastener
12" - 13" (35" Long)	1102009	1102008	1102007
14" - 19" (45" Long)	1004310	TAGTS1034751-NLA	1034744
20" - 24" (61" Long)	1034625	1034752-NLA See Note(s) A	1034745

## Miscellaneous Accessories

**Invacare does not service any parts associated with the TRRO/TRBKTS option including belts, belt attachment hardware, tie-down brackets, etc. Chairs with the Transport Ready or Transport Bracket option must be returned to Invacare for service.**

NOTE: A - 1034752 Airline Buckle Style is No Longer Available. Use Auto Style part number 1034625 as a replacement.  
Pelvic Belt is not available as a service part. Transport Ready chairs must be returned to Invacare for service with regard to specific Transport Ready Option parts.

## Miscellaneous Accessories

**Invacare does not service any parts associated with the TRRO/TRBKTS option including belts, belt attachment hardware, tie-down brackets, etc. Chairs with the Transport Ready or Transport Bracket option must be returned to Invacare for service.**

Part Number	Description	Package Quantity	Assembly Quantity
1027644-NLA	Calf Strap, 12" - 15" Wide Chairs (26" Long)	1	1
1032243	Calf Strap, 16" - 22" Wide Chairs (36" Long)	1	1
1163445-NLA	Calf Strap, Bodypoint 13" - 15" Wide Chairs (31-1/2" Long)	1	1
1163446	Calf Strap, Bodypoint 15" - 18" Wide Chairs (35-1/2" Long)	1	1
1163447	Calf Strap, Bodypoint 18" - 22" Wide Chairs (39-1/2" Long)	1	1
1032214	H-Strap, 12" - 15" Wide Chairs	1	1
1032215	H-Strap, 16" - 18" Wide Chairs	1	1
1032216	H-Strap, 19" - 22" Wide Chairs	1	1
1159758	Protector, Frame, Front, Blue	1	2
1159756	Protector, Frame, Back, Blue	1	1
1159384-NLA	Protector, Frame, Front, Silver	1	2
1159383-NLA	Protector, Frame, Back, Silver	1	1
1159759	Protector, Frame, Front, Red	1	2
1159757	Protector, Frame, Back, Red	1	1
1163422	Protector, Frame, Front, Black	1	2
1163421	Protector, Frame, Back, Black	1	1
1118578 See Note(s) A	Tilt Lock Out Lever Kit	1	2
1153995 See Note(s) B	Tilt Lock Out Collar Kit	1	2
See Positioning Strap Chart Category	Positioning Strap, Seat/Chest	1	See Positioning Strap Chart Category
See Stroller Handle, Non-Adjustable Category	Stroller Handle, Non-Adjustable	1	See Stroller Handle, Non-Adjustable Category

## Miscellaneous Accessories

**Invacare does not service any parts associated with the TRRO/TRBKTS option including belts, belt attachment hardware, tie-down brackets, etc. Chairs with the Transport Ready or Transport Bracket option must be returned to Invacare for service.**

Part Number	Description	Package Quantity	Assembly Quantity
See Stroller Handle, Adjustable Angle Category	Stroller Handle, Adjustable Angle	1	See Stroller Handle, Adjustable Angle Category
See Rigidizing Bar Category	Rigidizing Bar, Adjustable Width	1	See Rigidizing Bar Category
See Anti-Tipper Category	Anti-Tipper	1	See Anti-Tipper Category
See Bodypoint Side-Release Ankle Huggers & Hardware Category	Bodypoint Side-Release Ankle Huggers & Hardware	2	See Bodypoint Side-Release Ankle Huggers & Hardware Category
See Bodypoint Anterior Trunk Support Harness and Hardware Category	Bodypoint Anterior Trunk Support Harness and Hardware	1	See Bodypoint Anterior Trunk Support Harness and Hardware Category
See Back Pack with Mounting Straps Category	Back Pack with Mounting Straps	1	See Back Pack with Mounting Straps Category
See Oxygen Holder Category	Oxygen Holder	1	See Oxygen Holder Category
See I.V. Telescoping Rod Category	I.V. Telescoping Rod	1	See I.V. Telescoping Rod Category

NOTE: A - Includes Lever and attaching hardware. (Option RLL)  
 B - Includes Tilt Lock Out Collars and and Set Screw (Option TLC)

Pelvic Belt is not available as a service part. Transport Ready chairs must be returned to Invacare for service with regard to specific TRRO parts.